## **REMOWELL2 CVR**

*Cardiotomy/Venous Reservoir  
for extracorporeal circulation in cardio surgery*

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## **A.L.ONE AF PLUS**

*Oxygenator integrated with heat exchanger  
for extracorporeal circulation in cardio surgery*

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
**MODE D'EMPLOI**

---

**Figure 1 - REMOWELL2 CVR  
Cardiotomy/Venous Reservoir  
(flow direction indicated by bold  
arrows)**

1	Venous blood inlet port (VENOUS RETURN) 360° rotatable
2	2 female luer lock ports on venous blood inlet
3	Temperature probe port on venous blood inlet
4	Venous blood outlet port
5	Intracavitary cardiomy turret 360° rotatable
6	2 x 1/4" (6.35 mm) filtered ports on intracavitary cardiomy
7	1 x 3/8" (9.53 mm) filtered port on intracavitary cardiomy
8	1 x female luer lock filtered port on intracavitary cardiomy
9	1 x female filtered POS lock (ISO 8838) reducible to female luer lock on intracavitary cardiomy
10	Extracavitary cardiomy turret 360° rotatable
11	2 x 1/4" (6.35 mm) filtered ports on extracavitary cardiomy
12	1 x 3/8" (9.53 mm) filtered port on extracavitary cardiomy
13	1 x female luer lock filtered port on extracavitary cardiomy (for use as prime port)
14	1 x female filtered POS lock (ISO 8838) reducible to female luer lock on intracavitary cardiomy
15	Selector for extracavitary cardiomy opening
16	Cell saver port
17	1 x unfiltered female luer lock port
18	VENT/VACUUM port (VAVD technique)
19	Pressure relief valve
20	Sampling manifold
21	Reservoir top
22	Reservoir housing
23	Recess for holder
24	Holder
25	Knob
26	Holder pin

Figure 2 – Venous reservoir volume graduated scale



Figure 3 – Extracavitary cardiotomy volume graduate scale

## A. DESCRIPTION

REMOWELL2 CVR is a hard-shell cardiotomy/venous reservoir integrated with two cardiotomy filters, designed to allow venous drainage of the patient's blood, both through the hydrostatic load (height difference between the patient and the reservoir) and the vacuum-assisted venous drainage (VAVD) technique.

The device includes three compartments:

1. Venous reservoir, for venous blood collection, storage and filtration,
2. Intracavitary cardiotomy compartment, for intracavitary blood filtration,
3. Extracavitary cardiotomy compartment, for extracavitary blood filtration and lipid and leukocytes reduction.

The extracavitary cardiotomy compartment is separated from the venous reservoir through a partition wall which can be kept closed or opened rotating the selector on the reservoir top to the lock or unlock position respectively marked with  and  symbols. By turning the selector to the unlock position , the filtered/treated extracavitary blood flows into the venous reservoir. Depending on these operating conditions, the volume inside the venous reservoir varies and thus starting from the 850 ml volume the venous reservoir shows a volume double graduation scale, one yellow coloured marked with the symbol  indicating the volume in the closed selector condition (separated extracavitary cardiotomy), the other blue coloured, marked with the symbol  indicating the volume

in the open selector condition (venous reservoir containing the extracavitary blood) (Fig. 2).

The distinctive characteristic of REMOWELL2 CVR is the reduction of lipidic particles and leukocytes in the portion of extracavitary blood, coming from the surgical field and gradually collected in the pericardial sac following the surgical operations. Reduction of lipidic particles and leukocytes in extracavitary blood is achieved through:

1. **Retention by cascade filtration process:** sequential filtering media layers enable to partially hold back leukocytes and lipidic particles,
2. **Decantation and separation of the supernatant:** thanks to the sedimentation of the blood collected in the extracavitary cardiotomy reservoir and to a special skimming system, further lipidic particles collecting in the supernatant can be removed.

The values of the lipids and leukocytes removed are strongly affected by haemodilution of the blood aspirated from the surgical field and by the sedimentation time in the extracavitary cardiotomy reservoir.

The device is fitted with a pressure relief valve.

The reservoir's inner contact surfaces are coated with A.G.I.L.E. (Advanced Generation Inert Layer E.C.C.) system, based on Phosphorylcholine (PC), improving the device blood compatibility by reducing platelet adhesion on the coated surface.

The device is single-use, non-pyrogenic, supplied **STERILE** and individually packed. Sterilised by ethylene oxide.

## B. TECHNICAL FEATURES

Max. Blood Volume Capacity Venous Reservoir	4500 ml
Max. Blood Volume Capacity Extracavitary Cardiotomy	800 ml
Max Operating Volume	4500 ml
Min. Operating Volume	150 ml up to 5 l/min
	200 ml from 5 to 8 l/min
Blood Flow Range	Venous Flow: 0.5-8 l/min
	Cardiotomies: max. 2.5 l/min (each)
	Combined Venous and Cardiotomies flow: max. 8 l/min
Dynamic Priming (Hold up) Volume	
Venous Filter	250 ml @ 8 l/min
Volume scales tolerance	± 10%
<b>Filtration</b>	
Venous	Pore size 80 µm
	Efficiency ≥90% at 80 µm
Cardiotomies	Pore size 40 µm
	Efficiency ≥90% at 40 µm
<b>Lipidic Particles Filtration Extracavitary Cardiotomy Reservoir</b>	
Reduction Efficiency @ 60 min. <sup>1</sup>	Average 61%

	Min. 40%
	Max. 82%
<b>Leukocytes Filtration Extracavitary Cardiomy Reservoir</b>	
Reduction Efficiency @ 60 min. <sup>1</sup>	Average 36%
	Min. 18%
	Max. 54%
<b>Pressure relief valve</b>	
Positive pressure	5.3 mmHg
Negative pressure	-75.2 mmHg
<b>Max negative reservoir pressure</b>	
Negative pressure	-75.2 mmHg
<b>Ports</b>	
Venous blood inlet port	1/2" (12.7mm) – 360° rotatable
Luer lock at venous inlet	2 x female luer lock
Venous temperature probe port at venous inlet	YSI Series 400 fitting
Intracavitary cardiomy turret (360° rotatable)	2 x 1/4" (6.35 mm)
ports (mainly blood with limited amount of air)	1 x 3/8" (9.53 mm)
	1 x female luer lock
	1 x female POS lock (ISO 8838) reducible to female luer lock
	2 x 1/4" (6.35 mm)
Extracavitary cardiomy turret (360° rotatable)	1 x 3/8" (9.53 mm)
ports (mainly air with limited amount of blood)	1 x female luer lock
	1 x female POS lock (ISO 8838) reducible to female luer lock
	1 x 1/4" (6.35 mm)
Cell saver port	1 x 1/4" (6.35 mm)
VENT/VAACUM port	1 x 1/4" (6.35 mm)
Unfiltered luer lock	1 x female luer lock
Sampling manifold	Arterial blood sampling line (red handle stopcock)
	Venous blood sampling line (blue handle stopcock)
	Additional sampling/injection port (white handle stopcock)
	3/8" (9.53 mm)
Venous blood outlet port	
<b>Adapters</b>	
Venous blood inlet port adapter	1 x 1/2" (12.7mm) / 3/8" (9.53 mm)
Cardioplegia adapter	1 x POS lock male (ISO 8837) with 1/4" (6.35 mm) connector
	1 x 3/8" (9.53 mm) / 1/4" (6.35 mm)
Other adapters	1 x 1/4" (6.35 mm) / female luer lock
	1 x 3/8" (9.53 mm) / female luer lock
<b>Materials</b>	
Housing parts (reservoir and lid)	Polycarbonate (PC)
Filtering media (venous and cardiomy)	Polyester (PET), Polypropylene (PP)
Defoamer (venous and cardiomy)	Silicone treated Polyurethane (PU)
Seals	Thermoplastic Elastomer - Styrene-Butilene-Etilene-Stirene (SEBS)
	Polycarbonate (PC), Polypropylene (PP), Thermoplastic Elastomer - Styrene-Butilene-Etilene-Stirene (SEBS), Polyoxymethylene (POM)
Caps, adapters, selector	Phosphorylcholine (PC)
Biocompatible coating	Blood cell damage
<b>Technical data<sup>2</sup> available on request</b>	Air-handling capability
	Antifoam characteristics
	Break-through volume

<sup>1</sup> Data obtained in in-vitro testing.

<sup>2</sup> The data include the applicable tolerances.

### C. INTENDED PURPOSE

REMOWELL2 CVR hard-shell cardiotomy/venous reservoir is intended to be used in an extracorporeal perfusion circuit to collect, store and filter venous and cardiotomy suctioned blood during routine cardiopulmonary bypass procedures up to 6 hours in surgery.

The hard-shell cardiotomy/venous reservoir is also intended for use during vacuum-assisted venous drainage (VAVD) procedures.

REMOWELL2 CVR is furthermore intended to reduce the content of lipidic particles and leukocytes in the portion of extracavitary blood coming from the surgical field and gradually collected in the pericardial sac following the surgical operations.

REMOWELL2 CVR should not be used for more than 6 hours. Contact with blood for longer periods is not advised.

REMOWELL2 CVR should be used in combination with the medical devices listed in section O (Medical devices for use with REMOWELL2 CVR).

### PATIENT POLULATION AND INDICATIONS FOR USE

REMOWELL2 CVR is indicated for use on adult patients undergoing extracorporeal circulation during cardiopulmonary bypass procedures.

### D. CONTRAINDICATIONS

Do not use the device for any purpose other than indicated.

Do not use if air leaks are observed during priming or ECC, this may result in air embolism to the patient and/or fluid/blood loss.

### E. RISKS AND SIDE EFFECTS

The possible risks and side effects of the REMOWELL2 CVR include:

- haemolysis,
- adverse blood / tissues reaction,
- coagulation/complement activation, platelet dysfunction,
- thrombosis,
- gaseous embolism,
- thromboembolism,
- blood loss,
- haemodilution,
- infection.

#### WARNING:

- **Any serious incident occurring in relation to the device should be reported to the manufacturer and the competent authority of the Member**

**State in which the user and/or patient is established.**

### F. SAFETY INFORMATION

Information to draw the user's attention to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text as follows:

#### WARNING

Indicates serious adverse reactions and potential safety hazards for the practitioner and/or the patient, which may occur during proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

#### CAUTIONS

Indicates that the user must take every possibly precaution for safe and effective use of the device.

The following is general safety information for the operator preparing to use the device. Specific safety information is also given in the instructions for use where this information is relevant for correct operation.

#### WARNINGS

-  **Carefully read the instructions before use.**
- **Store in a dry place between 0°C (32°F) and 40°C (104°F).**
- **FRAGILE, handle with care.**
- **Single-use device. Do not reuse. Do not clean or resterilize the device for subsequent reuse.**
- **Improper reuse may cause transfusion reactions, cross-infection and alteration of device performance.**
- **For use by professionally trained personnel only.**
- **The device must be used in accordance with these instructions. EUROSETS cannot be held responsible for damage deriving from inexperience or improper use.**
- **Carefully check the device seal during priming and ECC. Any leakage may give rise to loss of sterility, blood leakage, gaseous embolism. If a leakage during priming or ECC is noticed, replace the defective device following good perfusion practices.**
- **The connection lines must be correctly connected in order to prevent any potential tube constriction or occlusion which may lead to reduced blood and/or gas flow.**

- To safely apply vacuum to the device, carefully read section L of these instructions for use.
- The special pressure relief valve on the cardiotomy/venous reservoir allows to discharge positive pressures higher than 5.3 mmHg (0.7 kPa) and negative pressures lower than -75.2 mmHg (-10 kPa).
- After use, dispose of the device in accordance with the regulations in force on biologically hazardous waste in the country of use.
- For further information and/or in case of a complaint, contact EUROSETS or an authorised local representative.

### CAUTION

- Follow an anticoagulation protocol and monitor anticoagulation before, during and after the bypass.

The benefit of using an extracorporeal support must be weighed against the risk of systemic anticoagulation and it must be prescribed by a physician. Maintain adequate anticoagulation before and during the procedure.

- Federal law (USA) restricts this device to sale by or on the order of a physician.

## G.SETUP

### 1. POSITIONING THE HOLDER

#### 1.1. CVR HOLDER

Position the Eurosets CVR holder on the pump stand and secure it by means of the knob found at the end of the carrying arm.

#### 1.2. OXYGENATING SYSTEM (OXYGENATOR + CVR) HOLDER

Position the Eurosets oxygenating system (oxygenator + CVR) holder on the pump stand and secure it by means of the knob found at the end of the carrying arm.

Verify that the holder is in the right position (see Adult/Paediatric Oxygenators Holder Instructions for Use for two different height holder positions).

### 2. SECURING DEVICE ON THE HOLDER

#### WARNINGS

- Check the expiry date on the label. Do not use the device after this date.
- The device has been sterilized by ethylene oxide and is non-pyrogenic. Its sterility and integrity are guaranteed only if the package is intact.
- Carefully inspect the device before use. Do not use it if the package is wet, open or

damaged, if the device is visibly damaged or if the protective caps are not in place.

- The device must be used immediately after opening the sterile package.
- The device must be handled aseptically.
- Do not use solvents such as alcohol, ether, acetone, etc., as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This would cause damage which may compromise the integrity and proper functioning of the device.
- Avoid that the connectors accidentally come into contact with lubricating agents which may encourage accidental disconnection.

#### 2.1. REMOWELL2 CVR

Secure REMOWELL2 CVR on the holder.

Check that REMOWELL2 CVR fits perfectly and fully into place with the holder support inserts in the two housing holes.

#### 2.2. REMOWELL2 SYSTEM (OXYGENATOR + CVR)

Secure REMOWELL2 system (oxygenator + CVR) on the holder (see possible matches in Figure 4).

Check that REMOWELL2 system fits perfectly and fully into place with the holder support inserts in the two housing holes and check that the oxygenator is securely in place in the support jaws.

Check and, if necessary, tighten the connections of the lines connected to the oxygenator.

Only in this condition will it be possible to secure REMOWELL2 system by turning the fastening lever anticlockwise, until block is achieved.

### 3. CONNECTING THE CIRCUIT

#### WARNINGS

- All the connections downstream of the pump must be secured by means of ties.
- Remove the yellow label (labelled with "Remove before use") and the yellow cap from the VENT/VACUUM connector.
- Before use remove from the pressure relief valve the yellow tab (labelled with "Remove before use").
- The action of the pressure relief valve does not prevent from excessive negative pressure in the reservoir caused by sudden abnormal high vacuum level (i.e. vacuum regulator device failure).

- Do not for any reason occlude the external access hole of pressure relief valve.
- Do not connect priming line or aspiration lines to the cell saver port, since this connector is dedicated to filtered blood recovery only.

**VENOUS LINE:** Connect a 1/2" venous line to the reservoir venous inlet port marked "VENOUS RETURN". The venous inlet port can be rotated through 360° to find the most convenient position for the venous tube.

**PRIMING LINE:** Connect the priming line to the desired extracavitary turret port and turn the selector to the lock position  to close the extracavitary cardiomy reservoir.

#### WARNINGS

- In case of fast priming need, connect priming line to one of the intracavitary turret ports.

**ASPIRATION LINES:** After removing the protective caps from the cardiomy turrets filtered inlet ports on cardiomy/venous reservoir top (each cardiomy turret composed of two 1/4" inlets, one 3/8" inlet, one Luer-lock inlet and one luer POS-lock), connect the ends of the aspiration lines and turn the turrets orienting the filtered inlets towards the suction pumps. Blood suctioned from intracavitary field (e.g. ventricle aspirator or aortic root aspirator) shall be connected to intracavitary cardiomy turret. Connect the oxygenator purge line to the intracavitary cardiomy turret. Blood suctioned from extracavitary field (e.g. blood collected from pericardial sac) shall be connected to extracavitary cardiomy turret.

**PUMP LINE:** The pump segment must be fitted between the venous reservoir outlet port and the oxygenator venous inlet port, taking account of the pump rotation direction.

#### 4. SAMPLING MANIFOLD

The arterial/venous sampling manifold is positioned on the venous reservoir lid.

If you wish to take samples from a more remote position, remove the manifold from the fastening guide, fit it on the dedicated holder and unwind the tubing provided. In this position, arterial/venous sampling can be performed in a radius of about 1 meter.

After performing artero/venous sampling, Eurosets recommends maintaining the three

stopcocks of the manifold in the close position, thus to allow an artero/venous loop during ECC, permitting a representative sampling of the perfusion status.

#### CAUTION

The user should check for occlusions in tubing during set up.

#### 5. CONNECTING THE TEMPERATURE PROBE

The venous temperature probe port is positioned on the venous inlet port of the reservoir.

#### 6. CONNECTING THE VACUUM LINE

If performing vacuum-assisted drainage (VAVD), please read carefully section L of these instructions for use.

### H. PRIMING PROCEDURE

#### WARNINGS

- Do not apply the vacuum during this procedure.
- Do not use alcohol-based priming solution, as this would compromise proper functioning of the device.

#### 1. CLOSE THE VENOUS LINE

Clamp the venous line leading from the patient.

#### 2. PRIME THE VENOUS RESERVOIR AND CARDIOTOMIES

Close all the aspiration lines with clamps and secure with ties all aspiration lines connected to the turrets.

It is recommended to carry out the filling from the extra-cavity line inlet port (Fig. 1 point 10), to ensure:

- correct filling and wetting of the filter system and
- filling of the siphon system necessary for blood decantation.

Priming solution quantity shall consider achievement of intended haematocrit value, taking into account that:

- the 3/8" tube capacity is 72 ml/m;
- the 1/2" tube capacity is 127 ml/m.

Clamp the cardiomy/venous reservoir outlet.

#### 3. PRIME THE OXYGENATOR, FOLLOWING THE INSTRUCTIONS FOR USE OF THE OXYGENATOR IN USE

Remove the pump segment of the arterial pump head and position it at the same height as the cardiomy/venous reservoir.

Slowly open the clamp at the cardiomy/

venous reservoir outlet and fill the tube, gradually lowering it with respect to the level in the cardiotomy/venous reservoir, so that it is completely filled and all the air is sent to the oxygenator.

The oxygenator module is hence completely filled by gravity. Once oxygenator is completely filled, position the pump segment in the arterial pump and start recirculation.

#### WARNING

- Do not use pulse flow during priming.
- Sudden changes in flow rate during priming can pull air across the membrane into the blood pathway.
- Check the correct dosage of anticoagulant in the system before starting the bypass.
- EUROSETS recommends using the pump speed regulator to reduce or slowly stop the arterial flow.
- Do not use the pump on/off switch until the pump speed is zero.

## I. INITIATING BYPASS

### 1. OPENING THE ARTERIAL AND VENOUS LINES

First remove the clamp from the arterial line and then the clamp from the venous line. Start bypass with a blood flow appropriate for the patient's size.

Constantly check the blood level in the cardiotomy/venous reservoir.

If performing vacuum-assisted drainage (VAVD), please read carefully section L of these instructions for use.

### 2. CLOSING THE EXTRACAVITARY CARDIOTOMY

Prior to start aspiration of extracavitary blood, turn the selector to the lock position , to allow selective blood collection and filtration in extracavitary cardiotomy reservoir.

## J. DURING BYPASS

### 1. CHECKING THE VENOUS RETURN

If a higher venous return flow is necessary, lower both the oxygenator and the cardiotomy/venous reservoir with respect to the patient's position.

If performing vacuum-assisted drainage (VAVD), please read carefully section L of these instructions for use.

#### WARNINGS

- The ACT (Activated Coagulation Time) must always be greater than or equal to 480 seconds

in order to ensure adequate anticoagulation of the extracorporeal circuit.

- If administration of anticoagulant to the patient is necessary, use the luer connector of the central stopcock (white handle) on the sampling manifold.
- The cardiotomy/venous reservoir minimum operating volume is 200 ml at 8 l/min and 150 ml at 5 l/min. Minimum operating volume is yellow highlighted on the reservoir volume graduated scale. To allow sufficient reaction time in case of sudden venous flow stop, preventing gaseous emboli passing to the patient, it is recommended to maintain the operating level at the double of the minimum operating volume.
- It is recommended to monitor the cardiotomy/venous reservoir volume using a reservoir level sensor.

## 2. CARDIOTOMY/VENOUS RESERVOIR MANAGEMENT

If performing vacuum-assisted drainage (VAVD), please read carefully section L of these instructions for use.

If you need to add blood mass compensating fluids to the venous reservoir, use the 1/4" or the luer-lock or the Pos-lock ports located on the intracavitary cardiotomy turret.

For the correct reading of the volume contained inside the venous reservoir in the conditions with open or closed selector, refer to indications reported in section A. Description.

The effectiveness of lipidic particles and leukocytes reduction is closely related to the blood dwell time in the extracavitary cardiotomy reservoir which ranges approximately from 30 to 60 minutes.

In the case of particularly short CPB or if needing blood mass volume, the extracavitary cardiotomy reservoir selector can be turned to unlock position .

In this condition, reduction of the lipids and leucocytes will be less effective (reduced filtration time for leukocytes and decreased sedimentation effect for lipids).

The person responsible for conducting CPB needs to evaluate when using the extracavitary blood before the required dwell time elapsed, when in need of blood mass volume.

The blood collected in the extracavitary cardiotomy reservoir can be:

- Recovered in the venous reservoir by turning the selector to the unlock position .

- Treated separately recovering it from the cell saver port.

With the selector turned to the unlock position , the blood coming from the extracavitary cardiotomy reservoir mixes with that present in the venous reservoir. The skimming system will prevent passage of the surnatant portion (approximate amount of 50 ml).

#### WARNING

- **If you wish to change from vacuum-assisted drainage to conventional venous drainage please read carefully section L of these instructions for use.**

### K. TERMINATING BYPASS

This must be carried out depending on the condition of each individual patient.

To recover extracavitary cardiotomy blood:

- 1) remove the cap from the cell saver port and connect port to a cell saver machine for blood components separation through a 1/4" line;
- 2) process the blood collected with the cell saver and reinfuse it to the patient.

#### CAUTION

If extracorporeal circulation has to be restarted, maintain a minimum blood flow inside REMOWELL2 CVR venous reservoir.

### L. BLOOD RECOVERY AFTER BYPASS

- 1) Recover all the blood contained in the venous line into the cardiotomy/venous reservoir as soon as the surgeon has removed the cannula from the patient's vena cava.
- 2) Perfuse through the aortic cannula as required by the patient's condition, slowly decreasing the level in the cardiotomy/venous reservoir.
- 3) If desired, blood in the oxygenator may also be returned by adding clear prime to the venous reservoir when the blood in the venous reservoir has reached the minimum volume. Pump the prime slowly through the oxygenator ensuring that the venous reservoir never be emptied.
- 4) When the cardiotomy/venous reservoir is nearly empty, stop the arterial pump and clamp the arterial line.

### M. USE OF THE VACUUM-ASSISTED VENOUS DRAINAGE TECHNIQUE (VAVD)

The device can be used with vacuum-assisted drainage as it is fitted with a pressure relief valve.

To apply safely vacuum to the device, carefully read the following instructions for use.


#### 1. CONNECTING THE VACUUM LINE

If performing venous drainage using the vacuum-assisted technique (VAVD), the vacuum line must be connected to the VENT/VACUUM port. Ensure that label and yellow protection cap have been removed from VENT/VACUUM port.

Connect the Eurosets VAVD connection kit (figure 2) to it.

Connect the other end of the line to a vacuum regulator equipped with a condensation trap.

Figure 2 – Connection line to the aspiration source for VAVD

#### 2. DURING PRIMING, RECIRCULATION AND BYPASS

##### WARNINGS

- Do not apply the vacuum during priming and recirculation procedures.
- It is recommended not to exceed -75.2 mmHg [-10 KPa] when applying negative pressure to the cardiotomy/venous reservoir.
- Periodically check functioning of the vacuum regulating device and the degree of vacuum.
- During vacuum-assisted drainage, check that all the connections are closed and apply the vacuum to the cardiotomy/venous reservoir.
- The vacuum applied to the cardiotomy/venous reservoir must be adjusted in order to obtain adequate venous return.
- If you wish to change from vacuum-assisted drainage to conventional venous drainage, suspend the vacuum application by detaching the vacuum line so that the venous reservoir is brought back to atmospheric pressure.
- If extracavitary cardiotomy blood needs to be sent to the cell saver machine, remove the vacuum and allow return to atmospheric pressure before connecting the cell saver port to the cell saver machine line. Stop the

drainage to the cell saver machine before restarting the VAVD.

#### CAUTION

During the vacuum-assisted venous drainage procedure pay attention to involuntary opening of the Luer-locks or connectors positioned on the turrets, as this causes a reduction of the vacuum in the reservoir; open the connectors only for the time necessary for the operation required.

### 3. LOW FLOW RECIRCULATION

In all cases that are expected a circulatory arrest and the opening of a shunt between the arterial and venous circuit, to maintain blood on the move, it is necessary to stop the vacuum application.

If is necessary to interrupt vacuum application, remove one of the leakproof caps on the cardiotomy/venous reservoir and remove clamp from the collecting line to the suction source for the VAVD technique.

Resume the arterial pump flow to the value previously set. In case of vacuum-assisted drainage, re-insert leakproof cap on the cardiotomy/venous reservoir and re insert the clamp on the collecting line to the suction source. Gradually increase the negative pressure to reach an adequate blood drainage.

### N. USING THE CARDIOTOMY/VENOUS RESERVOIR FOR POST-OPERATIVE AUTOTRANSFUSION

To use the cardiotomy/venous reservoir for post-operative autotransfusion, proceed as follows:

- 1) Disconnect the purge lines.
- 2) Disconnect the arterial sampling line from the sampling manifold.
- 3) Separate the venous reservoir from the oxygenating module by pressing the release button and simultaneously turning the blue ring counterclockwise (shown by the open padlock).
- 4) Close all the ports of the cardiotomy/venous reservoir used during bypass, including the VENT/VACUUM port.

#### WARNING

- **When using the cardiotomy/venous reservoir for post-operative autotransfusion, pay attention when administering protamine because excessive dosage may cause thrombus formation in the cardiotomy/venous reservoir.**

## 0. DEVICE REPLACEMENT

### 1. REMOWELL2 CVR

A spare device must always be available during perfusion. After 6 hours of use with blood or if the perfusionist is of the opinion that the conditions are such that the safety of the patient may be compromised (insufficient performance, leaks, abnormal blood parameters, etc.), replace the cardiotomy/venous reservoir following the instructions below.

#### CAUTION

Use a sterile technique throughout the replacement procedure.

- 1) Suspend the application of vacuum to the venous reservoir (when vacuum-assisted drainage has been used).
- 2) Turn the main pump and the suction pump off.
- 3) Close the venous return and the blood outlet lines, using two clamps (5 centimetres apart).
- 4) Cut the venous return and the blood outlet lines in the section between the two clamps, leaving a sufficient length of tubing to allow reconnection.
- 5) Remove all the sampling lines, the suction lines, and in any case all the connection to the REMOWELL2 cardiotomy/venous reservoir.
- 6) Remove REMOWELL2 CVR from the holder.
- 7) Place a new REMOWELL2 CVR on the holder. Connect all the lines (inlet, outlet, other connection lines, vacuum line, sampling lines, etc.).
- 8) Verify all connections and secure with ties.
- 9) Prime New REMOWELL2 CVR as described in Section G Prime Procedure.
- 10) Remove the clamps from the venous and the arterial lines, close the oxygenator purge line and restart the bypass.
- 11) The residual blood in the replaced cardiotomy/venous reservoir may be recovered by connecting the outlet to the 3/8" inlet on the turret of the new cardiotomy/venous reservoir.

### 2. REMOWELL2 SYSTEM (OXYGENATOR + CVR)

A spare device must always be available during perfusion.

After 6 hours of use with blood or if the perfusionist is of the opinion that the conditions are such that the safety of the patient may be compromised (insufficient oxygenator performance, leaks, abnormal blood parameters, etc.), replace the

oxygenator following the instructions below.

### CAUTION

Use a sterile technique throughout the replacement procedure.

- 1) Suspend the application of vacuum to the venous reservoir (when vacuum-assisted drainage has been used).
- 2) Turn off the gas flow.
- 3) Close the venous return line using two clamps (5 centimetres apart).
- 4) Stop the arterial pump and close the arterial line using two clamps (5 centimetres apart) near the oxygenator.
- 5) Turn off the thermocirculator, clamp and remove the water lines.
- 6) Disconnect the gas line, all the monitoring and sampling lines, and the aspiration lines.
- 7) Cut the venous return line and the arterial line at the point between the two clamps, leaving a sufficient length of tubing to allow reconnection.
- 8) Remove REMOWELL2 system from the holder and remove the pump segment from the arterial pump.
- 9) Place a new REMOWELL2 system on the holder. Connect all the lines (venous line to the venous reservoir, arterial and gas lines to the oxygenator, pump segment from the CVR to the oxygenator, and possible vacuum line to the CVR.

### WARNING

- **In this phase, keep the venous and arterial lines clamped.**

- 10) Connect the water lines, turn on the thermocirculator and check the integrity of the new REMOWELL2 system.
- 11) Fill the extracavitary cardiotomy of the new REMOWELL2 system with priming solution through one of the turret ports (1/4" or 3/8").
- 12) Rotate the selector in the open position and wait some seconds for the filling of the venous reservoir.
- 13) Prime the new REMOWELL2 system and evacuate the micro air bubbles as described in the priming and recirculation procedure.
- 14) Check all the connections and secure them with ties.
- 15) Remove the clamps from the venous and the

arterial lines, close the purge line and restart the bypass.

- 16) Turn on the gas flow and adjust gas flow rate as required.
- 17) The residual blood in the replaced CVR may be recovered by connecting the outlet to the 3/8" inlet on the intracavitary cardiotomy turret of the new venous reservoir.
- 18) The blood contained in the oxygenator and in the heat exchanger may be recovered by connecting the arterial line outlet to the 3/8" inlet on the intracavitary cardiotomy turret of the new venous reservoir.

### P. MEDICAL DEVICES FOR USE WITH REMOWELL2 CVR

The device can be integrated into any perfusion system.

At present EUROSETS is not aware of any contraindication to use of the device with occlusive pumps or centrifugal pumps. The use of other types of pumps must be agreed with EUROSETS.

The circuit connections must be made using tubes of a diameter compatible with the dimensions of the ports on the device (3/8").

When using the vacuum-assisted drainage (VAVD), please read carefully section L of these instructions for use.

Eurosets recommends the use of YSI 400 series temperature probes.

### Q. LIMITED WARRANTY

These warranty conditions are in addition to any statutory rights of the purchaser pursuant to applicable law.

EUROSETS warrants that all reasonable care has been taken in the manufacture of this medical device as required by the nature of the device and the use for which the device is intended.

EUROSETS warrants that the medical device is capable of functioning as indicated in these instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

Nonetheless, EUROSETS cannot guarantee that incorrect diagnosis or therapy and/or the particular physical and biological characteristics of an individual patient do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected. EUROSETS, whilst emphasizing the need to adhere

strictly to the instructions for use and to adopt all the necessary precautions for proper use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from improper use of this device.

EUROSETS undertakes to replace the medical device in the event that it is defective at the time of placing on the market or, if shipped by EUROSETS, at the time of delivery to the end user, unless such defect is attributable to the purchaser.

The above replaces all other warranties, explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No EUROSETS representative, dealer, distributor or intermediary or any other industrial or commercial organisation is authorised to make any claim or provide additional warranties other than those expressly stated herein.


EUROSETS declines all responsibility for any modifications to the Warranty Conditions and the information/instructions for use expressly stated herein.

The purchaser undertakes to comply with the terms of this warranty and agrees, in the event of a dispute or litigation with EUROSETS, not to make claims based on alleged or proven changes or alterations made to the warranty by anyone contrary to and/or in addition to what has been agreed upon herein.

The relationship between the parties to the contract (even if not stipulated in writing) for which this warranty is granted, as well as any dispute related thereto or connected therewith, as well as any relationship or dispute regarding this warranty, its interpretation and execution, nothing excluded and/or reserved, are exclusively governed by Italian legislation and jurisdiction.


Any disputes shall exclusively be settled by the Court of Modena (Italy).

Figure 1 – A.L.ONE AF Plus oxygenator


- |    |  |
|----|--|
| 1  | Arterial blood sampling port, female luer lock, with one-way valve |
| 2  | Arterial blood outlet, 3/8" (9.53 mm)                              |
| 3  | Arterial temperature probe port                                    |
| 4  | Gas outlet, 3/8" (9.53 mm)   |
| 5  | Venous blood inlet, 3/8" (9.53 mm)                                 |
| 6  | Air purge line, 1/4" (6.35 mm)                                     |
| 7  | Cardioplegia port, 1/4" (6.35 mm)                                  |
| 8  | Water Hansen fitting (inlet/outlet)                                |
| 9  | Gas inlet, 1/4" (6.35 mm)  |
| 10 | Port for arterial monitoring, female luer lock                     |

Figure 2 – A.L.ONE AF Plus coupled with Eurosets CVR, assembled on holder

- |   |   |
|---|---|
| 1 | Console pole clamp                            |
| 2 | Console pole clamp knob                       |
| 3 | Height and rotation arm                       |
| 4 | Holder knobs (height and rotation adjustment) |
| 5 | Holder pole                                   |
| 6 | Cardiotomy/Venous Reservoir arm pin           |
| 7 | Cardiotomy/Venous Reservoir arm               |
| 8 | Oxygenator arm                                |
| 9 | Oxygenator locking lever                      |


**Figure 3 – Example of circuit system including A.L.ONE AF Plus oxygenator and Eurosets CVR**

- |    |                        |
|----|------------------------|
| 01 | VENOUS LINE            |
| 02 | PUMP LINE              |
| 03 | BLOOD LINE             |
| 04 | ARTERIAL LINE          |
| 05 | GAS MIXER              |
| 06 | GAS INLET LINE         |
| 07 | TEMPERATURE MONITOR    |
| 08 | WATER LINES (IN/OUT)   |
| 09 | HEATER-COOLER UNIT     |
| 10 | ARTERIAL PRESSURE LINE |
| 11 | PRESSURE MONITOR       |

Figure 5 – Connection and disconnection of A.L.One AF Plus oxygenator from CVR

## DESCRIPTION

A.L.ONE AF PLUS is a microporous hollow-fibre membrane oxygenator consisting of a gas exchange module with an integrated heat exchanger and an integrated 38µm arterial filter that ensures arterial blood filtration with removal of microaggregates and microemboli. A.L.One AF PLUS is equipped with an upper "valved" luer lock for arterial blood sampling, and a lower luer lock for blood pressure

monitoring. A.L.ONE AF PLUS inner contact surfaces are coated with A.G.I.L.E. (Advanced Generation Inert Layer E.C.C.) system, based on Phosphorylcoline (PC), improving the device blood compatibility by reducing platelet adhesion on the coated surface. The device is single use, non-pyrogenic, supplied **STERILE** and individually packed. Sterilised by ethylene oxide.

## TECHNICAL FEATURES

### OXYGENATOR MODULE

Recommended blood flow range

Recommended gas flow range

Max blood pressure

Max water pressure

Membrane type

Membrane surface area

Heat exchanger surface area

Arterial filter pore size

Arterial filter surface area

Air-handling capability

Static priming volume

(oxygenator module + heat exchanger)

Residual blood volume

### Ports

Oxygenator venous inlet

Oxygenator arterial outlet

Oxygenator gas inlet

Oxygenator gas outlet

Water ports

Arterial temperature probe port

Arterial sampling port

Cardioplegia line port

Connection for blood arterial pressure monitoring

Air-purge line port

Performance Data

0.5 - 7 l/min

0.5 - 14 l/min

750 mmHg (100 kPa)

1500 mmHg (200 kPa)

Microporous polypropylene

1.65 m<sup>2</sup>

0.08 m<sup>2</sup>

80 µm

38 µm

200 cm<sup>2</sup>

230 cm<sup>2</sup>

≥ 98%

225 ml

126 ml

3/8" (9.53 mm)

3/8" (9.53 mm)

1/4" (6.35 mm)

3/8" (9.53 mm)

1/2" Hansen coupling

YSI Series 400 fitting


female luer lock with one-way valve

1/4" (6.35 mm)


Female Luer Lock (lower than blood out)

1/4" (6.35 mm)


## Oxygen transfer rate


**Carbon dioxide transfer rate**

**Blood side pressure drop**

**Heat exchanger performance factor R**

\*ISO 7199 condition

Bovine blood PVC02= 45±5 mmHg

Hb= 12±1 g/dL

SV02= 65±5%

B.E.= 0±5 mEq/L

Temp= 37±1°C

Qb/Qg=1

FiO2= 100%

## Materials

Housing  
 Microporous hollow fiber membrane  
 Filtering media  
 Potting  
 Heat exchanger  
 Seals

Caps and adapters

Biocompatible coating

**Technical data<sup>1</sup> available on request**

Polycarbonate (PC)  
 Polypropylene (PP)  
 Polyester (PET)  
 Polyurethane (PU)  
 Stainless steel  
 Thermoplastic Elastomer - Styrene-Butylene-Ethylene-Styrene (SEBS)  
 Polypropylene (PP), Thermoplastic Elastomer - Styrene-Butylene-Ethylene-Styrene (SEBS)  
 Phosphorylcholine (PC)

Gas pathway pressure drop at max. blood:gas flow rates:  
 7 l/min : 14 l/min  
 Blood cell damage  
 Particle release

## Accessory lines

Purge/recirculation line  
 Drawing line for blood cardioplegia

<sup>1</sup> The data include the applicable tolerances.

## INTENDED USE AND INDICATIONS FOR USE

A.L.ONE AF PLUS is intended for use in extracorporeal perfusion circuit during cardiopulmonary bypass in cardiac surgery to oxygenate and remove carbon dioxide from the blood and regulate the blood temperature.

A.L.ONE AF PLUS, being integrated with arterial filter, enables filtration of arterial blood with removal of microemboli and microaggregates larger than 38 µm.

A.L.ONE AF PLUS is equipped with an upper "valved" luer lock for arterial blood sampling, and a lower luer lock for blood pressure monitoring.

The blood to be treated should contain anticoagulant.

A.L.ONE AF PLUS is an ADULT oxygenator.

A.L.ONE AF PLUS should not be used for more than 6 hours. Contact with blood for longer periods is not advised.

A.L.ONE AF PLUS should be used in combination with the medical devices listed in section "MEDICAL DEVICES FOR USE WITH A.L.ONE AF PLUS".

## CONTRAINDICATIONS

Do not use the device for any purpose other than indicated.

Do not use if air leaks are observed during priming or operation, this may result in air embolism to the patient and/or fluid/blood loss.

## RISKS AND SIDE EFFECTS

The possible risks and side effects of the A.L.ONE

AF PLUS include:

- haemolysis,
- adverse blood / tissues reaction,
- coagulation/complement activation, platelet dysfunction,
- thrombosis,
- gaseous embolism,
- thromboembolism,
- blood loss,
- haemodilution,
- infection.

## SAFETY INFORMATION

Information to draw the user's attention to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text as follows:

### WARNING

Indicates serious adverse reactions and potential safety hazards for the practitioner and/or the patient, which may occur during proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

### CAUTION

Indicates that the user must take every possible precaution for safe and effective use of the device. The following is general safety information for the operator preparing to use the device. Specific safety information is also given in the instructions for

use where this information is relevant for correct operation.

## WARNINGS

-  Carefully read the instructions before use.
- Store in a dry place between 0°C (32°F) and 40°C (104°F).
- **FRAGILE, handle with care.**
- **Single-use device. Do not reuse. Do not clean or resterilize the device for subsequent reuse.**
- **Improper reuse may cause transfusion reactions, cross-infection and alteration of device performance.**
- For use by professionally trained personnel only.
- The device must be used in accordance with these instructions. EUROSETS cannot be held responsible for damage deriving from inexperience or improper use.
- A.L.One AF PLUS is provided with two **NON-INTERCHANGEABLE** luer locks: a superior “valved” luer for arterial blood sampling, and a lower luer lock for blood pressure monitoring.
- Carefully check the device seal during priming and operation. Any leakage may give rise to loss of sterility, blood leakage, gaseous embolism. If you notice leakage during priming or operation, replace the defective device following good perfusion practices.
- During use, a spare A.L.ONE AF PLUS oxygenator must always be available.
- The connection lines must be correctly connected in order to prevent any potential tube constriction or occlusion which may lead to reduced blood and/or water and/or gas flow.
- Before the cardiopulmonary bypass, administer to the patient an adequate dosage of anticoagulant, and adjust dosage through monitoring during and after the bypass.
- It is recommended not to expose the blood to a temperature lower than 18°C (64°F) and higher than 37°C (99°F).
- Ensure that the venous reservoir outlet is always located over the oxygenator.
- The extracorporeal circulation shall be carefully and continuously checked throughout the cardiopulmonary bypass procedure.
- After use, dispose of the device in accordance with the regulations in force on biologically hazardous waste in the country of use.
- For further information and/or in case of a

**complaint, contact EUROSETS or an authorised local representative.**

## SET UP

### 1. POSITIONING THE HOLDER

#### 1.1. OXYGENATOR MODULE HOLDER

Position the A.L.ONE AF PLUS holder on the pump stand and secure it by means of the knob found at the end of the carrying arm.

Verify that the holder is in the right position (see Instructions for use of Adult/Paediatric oxygenator module holder for the two different height holder positions).

#### 1.2. OXYGENATING SYSTEM HOLDER (OXYGENATOR + CVR)

Position the oxygenating system (oxygenator + CVR) holder on the pump stand and secure it by means of the knob found at the end of the carrying arm.

Verify that the holder is in the right position (see Instructions for use of Adult/Paediatric oxygenator module holder for the two different height holder positions).

### 2. SECURING ON THE HOLDER

#### WARNINGS

- Check the expiry date on the label. Do not use the device after this date.
- The device has been sterilized by ethylene oxide and is non-pyrogenic. Its sterility and integrity are guaranteed only if the package is intact.
- Carefully inspect the device before use. Do not use it if the package is wet, open or damaged, if the device is visibly damaged or if the protective caps are not in place.
- The device must be used immediately after opening the sterile package.
- The device must be handled aseptically.
- Do not use solvents such as alcohol, ether, acetone, etc., as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This would cause damage which may compromise the integrity and proper functioning of the device.
- Avoid that the ports accidentally come into contact with lubricating agents which may encourage accidental disconnection.

#### 2.1. A.L.ONE AF PLUS

Remove the device from the sterile packaging.

Secure A.L.ONE AF PLUS on the holder and check that it is securely in place in the support jaws.

Check and, if necessary, tighten the connections of the lines connected to the oxygenator.

## 2.2. OXYGENATING SYSTEM (OXYGENATOR + CVR)

Remove the device from the sterile packaging. Secure the oxygenating system (oxygenator + CVR) on the holder (see possible matches in Figure 4).

Check that holder support inserts fits perfectly and fully in the two CVR housing holes and check that the oxygenator is securely in place in the support jaws (where present).

Check and, if necessary, tighten the connections of the lines connected to the oxygenator.

Only in this condition it will be possible to secure the oxygenating system (oxygenator + CVR) by turning the fastening lever anticlockwise (where present), until block is achieved.

## 3. THERMOCIRCULATOR SET UP

Connect the water tubes to the A.L.ONE AF PLUS Hansen fittings (indifferently to the water inlet or outlet).

### WARNINGS

- **The water temperature at the heat exchanger inlet must not exceed 42°C (108°F).**
- **The water pressure in the heat exchanger must not exceed 1500 mmHg (200 kPa).**

## 4. CHECK THE HEAT EXCHANGER

### WARNING

- **Check the heat exchanger by recirculating water inside the heat exchanger for a few minutes. The integrity of the unit is guaranteed if there are absolutely no leaks from the water compartment.**

## 5. CIRCUIT CONNECTIONS

### WARNINGS:

- **All the connections downstream of the pump must be secured by means of ties.**

**ARTERIAL LINE:** Remove the red cap on the arterial outlet of the oxygenator marked "ARTERIAL OUTLET" and connect a 3/8" line.

**PUMP LINE:** The pump segment must be fitted between the venous reservoir outlet port and the oxygenator venous inlet port taking account of the pump rotation direction.

**CARDIOPLEGIA LINE:** If oxygenated blood is necessary for blood cardioplegia, connect the A.L.ONE AF PLUS cardioplegia blood line to the

cardioplegia circuit.

### OXYGENATOR MODULE PURGING LINE:

Connect the purging line terminal to the venous reservoir inlet.

**ARTERIAL SAMPLING LINE:** Remove the protective cap from the Luer connector positioned on the arterial outlet and connect to it the arterial sampling line.

### WARNINGS:

- **The female luer arterial sampling port, located over the arterial outlet port, pre-connected with sampling manifold, cannot be used as outlet pressure measurement port due to the presence of integrated one-way valve.**
- **It is possible to measure blood pressure by connecting to the luer lock located under the arterial outlet and identified with a label.**
- **Male luer fitting not supplied with EUROSETS products may damage the one-way valve placed inside the oxygenator arterial sampling luer port.**

## 6. CONNECTING THE TEMPERATURE PROBES

The arterial temperature probe port is positioned near the oxygenator module arterial outlet, while the venous temperature probe port is positioned on the venous return connector of the CVR.

## 7. CLOSING THE PURGE LINE

Close the clamp on the purge line.

## 8. CLOSING THE CARDIOPLEGIA LINE

Close the clamp on the cardioplegia line.

Connect the cardioplegia line to the cardioplegia kit if in use.

## 9. CONNECTING THE GAS INLET AND OUTLET LINES

Remove the green cap from the gas inlet port marked "GAS INLET" and connect the 1/4" gas line.

Ensure that the gas supply is from a suitable air/oxygen mixer such as Sechrist or a system with similar technical features.

The special configuration of the "GAS ESCAPE" port prevents potential occlusion of the gas outlet.

A capnograph can be attached to this port.

### WARNINGS:

- **The "GAS ESCAPE" system has been designed to avoid any possible risk of occluding the gas outlet, which would cause**

immediate passage of air into the blood compartment.

- The user should check for occlusions in tubing during set up.

## PRIMING AND RECIRCULATION PROCEDURE WARNING

- Do not use alcohol-based priming solution, as this would compromise proper functioning of the oxygenator module.

### 1. KEEP THE GAS FLOW OFF

### 2. KEEP THE OXYGENATOR PURGE LINE AND CARDIOPLEGIA LINE CLOSED

Check that the purge line clamp and cardioplegia line clamp are closed.

### 3. CLOSE VENOUS AND ARTERIAL LINES

Clamp the venous line leading from the patient. Clamp the arterial line a few centimetres away from the oxygenator arterial outlet port.

### 4. CHECK THE HEAT EXCHANGER

Verify again the integrity of the heat exchanger, paying particular attention to possible water leaks.

### 5. PRIME THE VENOUS RESERVOIR

Fill the CVR with sufficient liquid to ensure that the intended hematocrit is obtained. Clamp the CVR outlet.

### 6. CIRCUIT PRIMING WARNINGS

- The pressure in the blood compartment of the oxygenator module must not exceed 750 mmHg (100 KPa).
- Do not occlude the gas escape port to avoid that pressure achieved into the gas side compartment becomes higher than pressure into the blood site compartment.
- The pressure in the blood compartment must always be higher than that in the gas compartment. This is to prevent the formation of gas emboli in the blood compartment.
- Be sure to remove all the gas emboli before starting the ECC (extracorporeal circulation); gas emboli are dangerous for patient.
- It is always advisable to use a pre-bypass filter during priming.

Remove the pump segment of the arterial pump head and position it at the same height as the cardiotomy/venous reservoir.

Slowly open the clamp at the reservoir outlet and fill the tube, gradually lowering it with respect to the level in the cardiotomy/venous

reservoir, so that it is completely filled and all the air is sent to the oxygenator.

The oxygenator module is hence completely filled by gravity. When A.L.ONE AF PLUS is full, position the pump segment in the arterial pump.

### 7. OPEN THE VENOUS AND ARTERIAL LINES

Remove the clamps from the venous and arterial lines and increase the flow rate to 1000 ml/min.

### 8. OPEN THE PURGE LINE

Check that the purge line is connected to the turret of the CVR.

Increase the arterial pump speed to the maximum flow rate.

Open the purge line of the oxygenator module in order to evacuate the air towards the CVR.

Ensure the filling and the de-bubbling of the cardioplegia line.

### 9. PURGE THE AIR CONTAINED IN THE CIRCUIT

During this phase it is necessary to tap the entire circuit in order to facilitate the removal of microbubbles from the tube walls. After some minutes in which the flow is maintained at a high rate, all air will be evacuated.

### 10. PRIME THE SAMPLING MANIFOLD AND THE BLOOD PRESSURE MONITORING LINE

Prime the sampling lines (arterial/venous) until all the air has been evacuated. Turn the stopcock selectors of the sampling manifold so as to fill automatically the arterial and venous line.

Fill the blood pressure monitoring line ensuring that there is a sufficient "meniscus" to prevent air from entering during monitoring.

### 11. PRIME THE CARDIOPLEGIA LINE

Ensure the filling and the de-bubbling of the cardioplegia line.

### 12. CLOSE THE PURGE / RECIRCULATION LINE

After 3-5 minutes of circulation at high flow, all the residual air will have been evacuated and the purge/recirculation line can be closed.

### 13. CLOSE THE VENOUS AND ARTERIAL LINES

After having verified the absence of air inside the circuit, stop the pump and clamp the venous and arterial line.

### WARNINGS

- Do not use pulsatile flow during priming.
- Sudden changes in flow rate during priming can pull air across the membrane into the blood pathway.
- Check the correct dosage of anticoagulant in the system before starting the bypass.

- EUROSETS recommends using the pump speed regulator to reduce or slowly stop the arterial flow.
- Do not use the pump on/off switch until the pump speed is zero.
- Do not turn off the heater-cooler.
- If the cardioplegia line has been connected, make sure that all the air has been evacuated.
- Do not apply negative pressure to the coronary outlet. Negative pressures in the blood compartment might cause the formation of gaseous microemboli.

## INITIATING BYPASS

### 1. OPENING THE ARTERIAL AND VENOUS LINES

First remove the clamp from the arterial line and then the clamp from the venous line. Start bypass with a blood flow appropriate for the patient's size. Before starting the procedure, check that anticoagulation level (ACT) is appropriate.

Constantly check the blood level in the cardiotomy/venous reservoir.

### 2. CHECKING PROPER FUNCTIONING OF THE HEAT EXCHANGER

Check the temperature of the venous and arterial blood.

### 3. SELECTING THE APPROPRIATE GAS FLOW

After initiating cardiopulmonary bypass, check the content of gas in the blood with a blood gas analysis and make any necessary corrections.

The regulation of gas flow, blood flow and  $\text{FiO}_2$  % in normothermia must be performed according to the table below:

Blood flow l/min	$\text{FiO}_2$ %	Gas:Blood flow
<3.5	50	1:1
3.5 – 4.8	60-70	1:1
>4.8	70-80	1:1

#### WARNINGS

- Always open the gas flow after the blood flow. The gas blood flow ratio must never exceed 2:1.
- Do not occlude the gas escape port to avoid that pressure achieved into the gas side compartment becomes higher than pressure into the blood site compartment.
- The pressure in the blood compartment must always be higher than that in the gas compartment. This is to prevent the formation of gas emboli in the blood compartment.
- During the heating and cooling phases, pay attention to the rate of temperature increase

and decrease (gradient).

## 4. BLOOD GAS ANALYSIS

After a few minutes of bypass, measure the gas content in the blood. Depending on the values found, adjust as follows:

High $\text{pO}_2$	==>	decrease $\text{FiO}_2$ %
Low $\text{pO}_2$	==>	increase $\text{FiO}_2$ %
High $\text{pCO}_2$	==>	increase gas flow
Low $\text{pCO}_2$	==>	decrease gas flow

## DURING BYPASS

### 1. CHECKING THE VENOUS RETURN

If a higher venous return flow is necessary, lower both the oxygenator and the CVR respect to the patient's position.

With vacuum assisted drainage, the vacuum applied to the CVR must be adjusted in order to obtain adequate venous return.

#### WARNINGS

- The ACT (Activated Coagulation Time) must always be greater than or equal to 480 seconds in order to ensure adequate anticoagulation of the extracorporeal circuit.
- The CVR connected to the oxygenator shall be always placed over the oxygenator.
- If administration of anticoagulant to the patient is necessary, use the luer port of the central stopcock on the sampling manifold.
- To allow sufficient reaction time in case of sudden venous flow stop, preventing gaseous emboli passing to the patient, it is recommended to maintain an operating level inside the CVR double of the minimum operating volume declared.

### 2. ARTERIAL SAMPLING

For arterial blood sampling, connect the extension tubes of the sampling tap to the female Luer-lock connector (positioned on the arterial outlet). This Luer connector is fitted with a one-way valve for suction only.

#### WARNINGS

- Perform sampling only when the pump is running, otherwise the blood pressure in the compartment would drop and might lead to the formation of air bubbles.

### 3. LOW FLOW RECIRCULATION

(Hypothermia associated with circulatory arrest)

#### WARNINGS

- Suspend vacuum application throughout the

**circulatory arrest procedure.**

- **Bleed the vacuum line to bring the CVR back to atmospheric pressure by removing the vacuum connection.**

- 1) Reduce the gas flow rate to less than 500 ml/min,
- 2) Open the recirculation line of the arteriovenous loop.
- 3) Reduce the flow rate of the arterial pump to 500 ml/min and clamp the venous reservoir inlet line,
- 4) Clamp the oxygenator arterial line after the recirculation loop and recirculate at a maximum flow rate of 500 ml/min throughout the patient's circulatory arrest
- 5) To restart bypass after circulatory arrest, open the venous and arterial lines and slowly increase the blood flow,
- 6) Close the recirculation line used and adjust the gas flow.
- 7) Resume the arterial pump flow to the value previously setted.

#### **WARNINGS**

- **Do not use purging line during low flow recirculation because oxygenator module doesn't work properly.**

#### **4. CONTINUOUS AIR PURGING**

If air continuously arrives at the oxygenator, hold the clamp of the purge line partially open so as to eliminate the air (this procedure however entails a drop in flow rate to the patient, proportional to the clamp opening, to the arterial pump flow rate and to the arterial pressure).

#### **TERMINATING BYPASS**

This must be carried out depending on the condition of each individual patient. Act as follows:

- 1) Turn off the gas flow.
- 2) Turn the heater-cooler off.
- 3) Slowly reduce the arterial flow to zero, at the same time closing the venous line.
- 4) Open the recirculation line.
- 5) Occlude the arterial line after the recirculation zone.
- 6) Increase the arterial flow rate to 500 ml/min.

#### **WARNINGS**

- **If extracorporeal circulation has to be restarted, maintain a minimum blood flow inside A.L.ONE AF PLUS (maximum 1500 ml/min).**

- **During recirculation do not turn the thermocirculator off.**
- **Check that the cardioplegia circuit connected to the coronary outlet is properly clamped.**
- **We recommend that you disconnect the CVR from the oxygenator module when the same is still on the holder (lift up the lever and rotate counterclockwise the blue ring, indicated by the lock open symbol).**

#### **RECOVERING THE BLOOD VOLUME AFTER BYPASS**

- 1) Recover all the blood contained in the venous line into the CVR as soon as the surgeon has removed the cannula from the patient's vena cava.
- 2) Perfuse through the aortic cannula as required by the patient's condition, slowly decreasing the level in the CVR.
- 3) If desired, blood in the oxygenator may also be returned by adding clear prime to the venous reservoir when the blood in the reservoir as reached minimum value. Pump the prime slowly (minimum flow ~50 ml/min) through the oxygenator ensuring that the reservoir never be emptied.
- 4) When the CVR is nearly empty, stop the arterial pump and clamp the arterial line.

#### **DEVICE REPLACEMENT**

A spare device must always be available during perfusion.

After 6 hours of use with blood or if the perfusionist is of the opinion that the conditions are such that the safety of the patient may be compromised (insufficient oxygenator performance, leaks, abnormal blood parameters, etc.), replace the oxygenator following the instructions below.

#### **1. A.L.ONE AF PLUS**

##### **CAUTION**

**Use a sterile technique throughout the replacement procedure.**

- 1) Turn off the gas flow.
- 2) Turn the arterial pump off and close, by means of a double clamp, the arterial line (5 centimetres apart) placed next to the oxygenator.
- 3) Close, by means of a double clamp, the venous line (5 centimetres apart).
- 4) Turn the thermocirculator off, clamp and remove the water lines.

- 5) Disconnect the gas line, all monitoring and sampling lines.
- 6) Cut the inlet venous line of the oxygenator and the arterial lines in the section between the two clamps, leaving a sufficient length of tubing to allow re-connection.
- 7) Remove the A.L.ONE AF PLUS from the holder
- 8) Place a new A.L.ONE AF PLUS on the holder. Connect all lines (arterial and gas lines to the oxygenator, pump line from CVR to oxygenator, arterial monitoring line, cardioplegia line and sampling line).
- 7) Cut the venous return line and the arterial line at the point between the two clamps, leaving a sufficient length of tubing to allow reconnection.
- 8) Remove the oxygenating system (oxygenator + CVR) from the holder and remove the pump segment from the arterial pump.
- 9) Place a new oxygenating system (oxygenator + CVR) on the holder. Connect all the lines (venous line to the CVR, arterial and gas lines to the oxygenator, pump segment from the CVR to the oxygenator, arterial monitoring line and possible vacuum line to the CVR).

### WARNINGS

- **In this phase, keep the venous and arterial lines clamped.**

- 9) Open the water lines, turn the thermocirculator on and check the integrity of the heat exchanger of the new A.L.ONE AF PLUS.
- 10) Fill the new A.L.ONE AF PLUS with priming solution and evacuate the microbubbles, as described in the priming and recirculation procedure.
- 11) Verify all connections and secure with ties.
- 12) Remove clamps from the venous and the arterial line, close the purging/recirculation line and start the bypass again.
- 13) The blood contained in the removed oxygenator may be recovered by connecting the arterial line to one of the inlet connectors of the venous reservoir.

## 2. OXYGENATING SYSTEM (OXYGENATOR + CVR) CAUTION

Use a sterile technique throughout the replacement procedure.

- 1) Suspend the application of vacuum to the CVR (when vacuum-assisted drainage is in use).
- 2) Turn off the gas flow.
- 3) Close the venous return line using two clamps (5 centimetres apart).
- 4) Stop the arterial pump and close the arterial line using two clamps (5 centimetres apart) near the oxygenator.
- 5) Turn off the thermocirculator, clamp and remove the water lines.
- 6) Disconnect the gas line, all the monitoring and sampling lines, and the aspiration lines.

### WARNING

- **In this phase, keep the venous and arterial lines clamped.**

- 10) Connect the water lines, turn on the thermocirculator and check the integrity of the new oxygenating system (oxygenator + CVR).
- 11) Fill the CVR of the new oxygenating system (oxygenator + CVR) with priming solution through one of the turret ports (1/4" or 3/8").
- 12) Prime the new oxygenating system (oxygenator + CVR) and evacuate the micro air bubbles as described in the priming and recirculation procedure.
- 13) Check all the connections and secure them with ties.
- 14) Remove the clamps from the venous and the arterial lines, close the recirculation/purge line and restart the bypass.
- 15) Turn on the gas flow and adjust gas flow rate as required.
- 16) The residual blood in the replaced CVR may be recovered by connecting the outlet to the 3/8" inlet on the turret of the new CVR.
- 17) The blood contained in the replaced oxygenator may be recovered by connecting the arterial line outlet to the 3/8" inlet on the turret of the new CVR.

### 2.1. REPLACING ONLY THE OXYGENATOR MODULE IN THE OXYGENATING SYSTEM

- 1) Suspend the application of vacuum to the CVR (when vacuum-assisted drainage is in use).
- 2) Close the gas flow and disconnect the gas line.

- 3) Clamp the venous return line.
- 4) Stop the arterial pump and close the arterial line with two clamps (5 centimetres apart) near the oxygenator.
- 5) Close the oxygenator inlet line near the oxygenator module port with two clamps (5 centimetres apart).
- 6) Turn off the thermocirculator, clamp and remove the water lines.
- 7) Close the purge clamp and the cardioplegia line using a metal clamp where necessary.
- 8) Cut the lines at the point between the two clamps.
- 9) Make sure that the sampling manifold is closed and disconnect the arterial sampling line being careful not to contaminate the luer lock of the sampling manifold (e.g. by connecting it to a female luer lock positioned on the top of the reservoir).
- 10) Remove all the other monitoring lines.
- 11) Cut the oxygenator inlet line and the arterial line at the point between the two clamps, leaving a sufficient length of tubing to allow reconnection.
- 12) Turn the holder locking lever (if present).
- 13) Lift up the lever present on the CVR base and rotate counterclockwise the blue ring (indicated by the lock open symbol) to unlock, as per Figure 5.
- 14) Remove the oxygenator module.
- 15) Position a new oxygenator module on the holder (if present).
- 16) Lock the new oxygenator module to the CVR. The CVR base shall be inserted in the blue ring and rotated clockwise until the lever locks into the slot, as per Figure 5.
- 17) Connect the oxygenator inlet line and the arterial line to the oxygenator module.
- 18) Connect and turn on the thermocirculator.
- 19) If possible, connect the reservoir to the oxygenator module by turning the blue ring nut (shown with the padlock closed).
- 20) Connect the purge line to the turret of the CVR; reconnect the cardioplegia line.
- 21) Connect the gas line, the arterial sampling line and all the monitoring lines.
- 22) Remove the clamp from the venous line.
- 23) Open the recirculation line (filter purge or arteriovenous loop), fill the oxygenator module by means of the arterial pump, eliminate the air in the module by opening the purge line in sequence, and recirculate

at a maximum flow rate of 1500 ml/min.

- 24) Close the recirculation line remove the clamp from the arterial line and start the new bypass.
- 25) Close the purge line and reconnect the cardioplegia circuit if necessary.
- 26) Restore the vacuum in the CVR if using the vacuum-assisted drainage technique.

## **MEDICAL DEVICES FOR USE WITH A.L.ONE AF PLUS**

The use of A.L.ONE AF PLUS is bound by the use of a hard-shell venous reservoir (open system) or a collapsible venous reservoir (closed system) and by the use of a cardiotomy reservoir for the recovery of blood from the surgical field.

Furthermore, an arterial / venous sampling system must be used.

At present EUROSETS is not aware of any contraindication to use of the device with occlusive pumps or centrifugal pumps. The use of other types of pumps must be agreed with EUROSETS.

Use Sechrist or a system with similar technical features as air/oxygen mixer.

Eurosets recommends the use of YSI temperature probes.

The circuit connections must be made using tubes of a diameter compatible with the dimensions of the ports on the device (3/8", 1/4", 1/2").

There are no contraindications for use of any heating/cooling system (thermocirculator).

## **LIMITED WARRANTY**

These warranty conditions are in addition to any statutory rights of the purchaser pursuant to applicable law.

EUROSETS warrants that all reasonable care has been taken in the manufacture of this medical device as required by the nature of the device and the use for which the device is intended.

EUROSETS warrants that the medical device is capable of functioning as indicated in these instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

Nonetheless, EUROSETS cannot guarantee that incorrect diagnosis or therapy and/or the particular physical and biological characteristics of an individual patient do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

EUROSETS, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the necessary precautions for proper use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from improper use of this device.

EUROSETS undertakes to replace the medical device in the event that it is defective at the time of placing on the market or, if shipped by EUROSETS, at the time of delivery to the end user, unless such defect is attributable to the purchaser.

The above replaces all other warranties on sale and/or functionality.

No EUROSETS representative, dealer, distributor or intermediary or any other industrial or commercial organisation is authorised to make any claim or provide additional warranties other than those expressly stated herein.

EUROSETS declines all responsibility for any modifications to the Warranty Conditions and the information/instructions for use expressly stated herein.

The purchaser undertakes to comply with the terms of this warranty and agrees, in the event of a dispute or litigation with EUROSETS, not to make claims based on alleged or proven changes or alterations made to the warranty by anyone contrary to and/or in addition to what has been agreed upon herein.

The relationship between the parties to the contract (even if not stipulated in writing) for which this warranty is granted, as well as any dispute related thereto or connected therewith, as well as any relationship or dispute regarding this warranty, its interpretation and execution, nothing excluded and/or reserved, are exclusively governed by Italian legislation and jurisdiction.

Any disputes shall exclusively be settled by the Court of Modena (Italy).

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