

**Enhanced convenience and unparalleled safety.**



***THERASAFE***<sup>®</sup>

**Therapeutic Phlebotomy Solutions**

---

 **NIPRO**  
CANADA  
Live Longer. Live Better.

[www.nipro.ca](http://www.nipro.ca)

## Therapeutic Phlebotomy Solutions



The CM-TS17X1B and CM-S17X1BH both have 17G x 1" safety needle with "Back eye" to optimize flow and minimize the risk of aspiration to the vessel wall. [2]

The CM-TS18X1B and CM-S18X1BH both have 18G x 1" safety needle with "Back eye" to optimize flow and minimize the risk of aspiration to the vessel wall. [2]




Introducing the sampling and infusion connector, which is swabable and designed with a split septum.



Built in safety needle

New patient improvement initiative with our 18G Safetouch™ II safety needle. [1]

Our stainless steel phlebotomy bag stand (CM-PH1) can be easily cleaned after use.



## Ordering Information

Product Code	Description	Qty. per Box
CM-TS17X1B	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor & 17G x 1" Needle	50
CM-TS17X1B10	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor & 17G x 1" Needle	10
CM-TS17X1BH	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor, BCTH & 17G x 1" Needle	50
FS+172530BC	17G x 1" Safetouch™ II Fistula Needle w/Backeye, Clamp & Fixed Wing	50
CM-TS18X1B	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor & 18G x 1" Needle	50
CM-TS18X1B10	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor & 18G x 1" Needle	10
CM-TS18X1BH	TheraSafe® Phlebotomy Bag w/Needleless Y-Adaptor, BCTH & 18G x 1" Needle	50
FS+182530BC-G	18G x 1" Safetouch™ II Fistula Needle w/Backeye, Clamp & Fixed Wing, Gamma	50
CM-PH1	Phlebotomy Bag Holder	1

### REFERENCES

- [1] Canadian Vascular Access Association (CVAA), Canadian Vascular Access & Infusion Therapy Guidelines, Pembroke, Ont.: Pappin Communications, 2019.  
 [2] National Kidney Foundation, KDOQI Guidelines, 2006.

Safetouch™ is a trademark of Nipro Corporation.



NIPRO CANADA CORPORATION  
 199 St. David Street, Lindsay, Ontario, Canada K9V 5K7  
 Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518  
 www.nipro.ca | service@nipro.ca

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE

NIPRO CORPORATION  
 3-26, Senriokashinmachi, Settsu,  
 Osaka, 566-8510, Japan

© 2023 Nipro Canada. All Rights Reserved.  
 TheraSafe\_Rev\_04  
 2023-11