

# HighFlow catheters


## Consistant performance in a demanding field

HighFlow sets the bar for tunneled dialysis catheters. Fabricated from a special carbothane formulation, the HighFlow offers a number of important material characteristics that:

- allow for thin walls which maximize inner lumen diameters
- offer strength while resisting kinking
- result in material which is thermosensitive and compliant enough to optimize patient comfort
- resist deterioration from alcohol and non-alcohol based antiseptics

HighFlow also provides a number of additional features:

- our wider cuff allows for a larger surface area for tissue in-growth thus providing a more secure anchoring of the catheter
- french size, catheter body length and priming volumes appear on the catheter to facilitate quick reference
- available in a choice of a straight or pre-curved "IJ" design



The introducer tray includes several key features:

- an improved tunneler flexible enough to allow excellent "shape-ability" yet provide sufficient strength for tunneling
- a superior safety scalpel which locks in place securely

# HighFlow catheters

## ORDERING INFORMATION

ITEM	DESCRIPTION	SIZE & LENGTH	TIP TO CUFF	PRIMING VOLUME
CM-15HFKB	IJ HighFlow Dialysis Kit	14.5F x 24 cm	19 cm	1.9 cc Arterial / 1.9 cc Venous
CM-19HFKB	IJ HighFlow Dialysis Kit	14.5F x 28 cm	23 cm	2.0 cc Arterial / 2.0 cc Venous
CM-23HFKB	IJ HighFlow Dialysis Kit	14.5F x 32 cm	27 cm	2.2 cc Arterial / 2.2 cc Venous
CM-27HFKB	IJ HighFlow Dialysis Kit	14.5F x 36 cm	31 cm	2.4 cc Arterial / 2.4 cc Venous
CM-24HFK	HighFlow Dialysis Kit	14.5F x 24 cm	19 cm	1.9 cc Arterial / 1.9 cc Venous
CM-28HFK	HighFlow Dialysis Kit	14.5F x 28 cm	23 cm	2.0 cc Arterial / 2.0 cc Venous
CM-32HFK	HighFlow Dialysis Kit	14.5F x 32 cm	27 cm	2.2 cc Arterial / 2.2 cc Venous
CM-36HFK	HighFlow Dialysis Kit	14.5F x 36 cm	31 cm	2.4 cc Arterial / 2.4 cc Venous
CM-40HFK	HighFlow Dialysis Kit	14.5F x 40 cm	35 cm	2.5 cc Arterial / 2.5 cc Venous
CM-55HFK	HighFlow Dialysis Kit	14.5F x 55 cm	50 cm	3.2 cc Arterial / 3.2 cc Venous
CM-75HFK	HighFlow Dialysis Kit	14.5F x 75 cm	70 cm	3.6 cc Arterial / 3.6 cc Venous
CM-15HFBIRV	IJ HighFlow Interventional Radiology Kit	14.5F x 24 cm	19 cm	1.9 cc Arterial / 1.9 cc Venous
CM-19HFBIRV	IJ HighFlow Interventional Radiology Kit	14.5F x 28 cm	23 cm	2.0 cc Arterial / 2.0 cc Venous
CM-23HFBIRV	IJ HighFlow Interventional Radiology Kit	14.5F x 32 cm	27 cm	2.2 cc Arterial / 2.2 cc Venous
CM-27HFBIRV	IJ HighFlow Interventional Radiology Kit	14.5F x 36 cm	31 cm	2.4 cc Arterial / 2.4 cc Venous
CM-24HFIRV	HighFlow Interventional Radiology Kit	14.5F x 24 cm	19 cm	1.9 cc Arterial / 1.9 cc Venous
CM-28HFIRV	HighFlow Interventional Radiology Kit	14.5F x 28 cm	23 cm	2.0 cc Arterial / 2.0 cc Venous
CM-32HFIRV	HighFlow Interventional Radiology Kit	14.5F x 32 cm	27 cm	2.2 cc Arterial / 2.2 cc Venous
CM-36HFIRV	HighFlow Interventional Radiology Kit	14.5F x 36 cm	31 cm	2.4 cc Arterial / 2.4 cc Venous
CM-40HFIRV	HighFlow Interventional Radiology Kit	14.5F x 40 cm	35 cm	2.5 cc Arterial / 2.5 cc Venous
CM-55HFIRV	HighFlow Interventional Radiology Kit	14.5F x 55 cm	50 cm	3.2 cc Arterial / 3.2 cc Venous

To order catheters only (with dilators) – delete the “K” from the product number

To order a Maximum Barrier version of a catheter kit, add “MB” to the order code\* .

Unless otherwise indicated, products listed are latex-free, ETO sterilized and packaged 5 per box

\*The Maximum Barrier version of a kit includes a patient cover, surgical gown, surgical mask, head cover and non-latex surgical gloves in addition to standard components of a High Flow catheter kit.

### HighFlow catheter kits include:

- 1 – 14.5F HighFlow dialysis catheter
- 1 – 15F peel-away sheath & dilator
- 1 – 8” malleable tunneler with sheath
- 1 – 0.038” x 70cm with 3 cm J guidewire with advancer
- 1 – 18G x 7cm GWI introducer needle
- 1 – 5cc retractable safety syringe with 25G x 1 1/2” needle
- 1 – 10cc syringe
- 1 – 25G x 3 1/2” Quincke needle
- 4 – 4x4” gauze
- 2 – plug caps
- 1 – suture

Note: Kit is double sterile packaged.

### HighFlow catheter IR kits include:

- 1 – 14.5F HighFlow dialysis catheter
- 1 – 8” malleable tunneler with sheath
- 1 – 0.038” x 70cm with 3 cm J guidewire with advancer
- 1 – 15F peel-away introducer with valve
- 1 – 18G x 7cm GWI introducer needle



### CardioMed Supplies Inc.

199 Saint David Street,  
Lindsay, ON K9V 5K7 Canada  
Tel: 705.328.2518 Toll-free: 800.387.9757  
Fax: 705.328.9747  
Web: [www.cardiomed.com](http://www.cardiomed.com) Email: [mail@cardiomed.com](mailto:mail@cardiomed.com)