



PERSONAL PROTECTIVE EQUIPMENT



NIPRO
CANADA

Live Longer. Live Better.

www.nipro.ca

PERSONAL PROTECTIVE EQUIPMENT

Face Shield (CM-FP-01)

- Featuring an anti-fog and anti-static coating, Nipro Canada's Face Shield is the perfect choice for clinicians.
- Provides an excellent area of coverage, preventing liquid splashes or sprays from reaching the wearer's eyes.
- Its open-back design allows the wearer to easily hear without restriction, and the flexible plastic shield offers a high degree of clarity without the typical obstructions found in other face shields of this nature.
- A comfortable, elastic headband allows for a one-size-fits-all configuration, and the soft hypo-allergenic foam pad helps to absorb moisture. Compatible with glasses.

Surgical Face Mask ATSM Level 3 (FM-44EE)

- 4-Ply/Layer Design: Extra protection against liquid penetration.
- 4-Fold Design: 15% larger breathing space vs 3-fold.



Surgical Face Mask w/Eye Shield ASTM Level 3 (FMS-44EPU)

- PET eye shield: Anti-fog and anti-splash.
- Mask specification equivalent to FM-44EE.



Medical Face Mask Tests and Requirements

Test	ASTM F2100-19			EN 14683:2019
	Level 1	Level 2	Level 3	Type IIR
Bacterial Filtration Efficiency (%)	≥95	≥98	≥98	≥98
Differential Pressure (mmH ₂ O/cm ²) (Pa/cm ²)	<5.0 <49.0	<6.0 <58.8	<6.0 <58.8	<6.1 <60.0
Submicron particulate filtration efficiency at 0.1 micron (%)	≥95	≥98	≥98	Not Required
Splash Resistance/Synthetic Blood Resistance mmHg	80	120	160	> 120 (16.0kPa)
Flame Spread	Class 1	Class 1	Class 1	*

* SEE EUROPEAN MEDICAL DIRECTIVE (2007/47/EC, MDD 93/42/EEC)

Ordering Information

Product Code	Description	Quantity
CM-FP-01	Disposable Face Shield	10 per box, 200 per case
CM-FP-02	Procedural Ear-loop Mask	50 per box, 2000 per case
FM-44EE	ASTM Level III Face Mask Colour Blue	50 per box, 2000 per case
FMS-44EPU	ASTM Level III Face Mask w/Shield	50 per box, 400 per case
FM-44ET	ASTM Level III Face Mask w/Ties	50 per box, 2000 per case
FM-KG-01	KN95, Non-woven Facemask, Large	40 per box, 800 per case

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE



NIPRO CANADA CORPORATION
199 St. David Street, Lindsay, Ontario, Canada K9V 5K7
Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518
www.nipro.ca | service@nipro.ca

NIPRO CORPORATION
3-26, Senriokashinmachi, Settsu,
Osaka, 566-8510, Japan

© 2024 Nipro Canada. All Rights Reserved.
FaceShield_Rev02
2024-01