



PARASiL[®]

Paracentesis Drainage Kits

Patient Comfort, Safety, Value



NIPRO
CANADA


Live Longer. Live Better.

www.nipro.ca

Paracentesis Drainage Kits

ParaSil® Paracentesis Drainage Tray (CM-002PCK)

Designed for patients who require ascites drainage, Nipro Canada's ParaSil® Paracentesis Drainage Tray contains all the sterilized components necessary for safe insertion of the drainage catheter.



Benefits Offered

- The catheter's soft and flexible silicone composition conforms to the peritoneal space and minimizes patient discomfort.
- The catheter's radiopaque stripe facilitates its placement using medical imaging.
- The catheter's multiple side holes improve drainage flow.
- Includes needleless swabable connector to maintain a closed system. Compatible with the standard luer connector offered on Nipro Canada's SmartTap® 2L and 4L drainage bags.

Sterilized ParaSil® Paracentesis Drainage Tray Includes:


- Silicone Drainage Catheter - 42cm, Single Cuff
- Barbed Female Luer Adapter
- Needleless Swabable Connector
- Syringe w/M.L.L
- Needle Block
- Introducer Needle
- "J" Guidewire w/Advancer
- Peel-Away Sheath Dilator
- Tunneler with Tapered Sleeve
- Side Pinch Clamp



Paracentesis Drainage Kits

ParaSil[®] Paracentesis Drainage Kit (CM-002PDK) and Procedure Kit (CM-002PBK)

Designed to facilitate the drainage of continuous malignant ascites. Suitable for use in a home, hospital, hospice, or long-term care facility.



CM-002PDK shown

Benefits Offered

- Comes with sterilized fenestrated drape, pleated towel, and transparent dressing for establishing safe, clean treatment.
- Sterilized SmartTap[®] Drainage Bag (CM-002PDK – Parasil[®] Drainage Kit only).
 - ◆ Suitable volume for treatment, 2L capacity.
 - ◆ Comes with graduation lines to visualize the amount of pathological fluid removed.
 - ◆ Tubing line comes with bed sheet clip, pinch clamp, and male luer connector.

Sterilized ParaSil[®] Paracentesis Drainage Kit and Procedure Kit Include:

- Blue Underpad
- Pleated White Towel
- Fenestrated Drape
- Gauze Pads
- Syringe w/M.L.L
- Adhesive Transparent Dressing
- 2L SmartTap[®] Drainage Bag (CM-002PDK only)
- Tube Occluding Forceps



Paracentesis Drainage Kits

Ordering Information

Product Code	Description	Quantity Per Box
CM-0020	SmartTap [®] Paracentesis Drainage Bag, 4L	25
CM-0022	SmartTap [®] Paracentesis Drainage Bag, 2L	25
CM-002PCK	Parasil [®] Drainage Tray w/Peritoneal Catheter	5
CM-002PDK	ParaSil [®] Drainage Kit	8
CM-002PBK	Parasil [®] Procedure Kit	8



NIPRO
CANADA

Live Longer. Live Better.

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE



NIPRO CANADA CORPORATION

199 St. David Street, Lindsay, Ontario, Canada K9V 5K7
Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518
www.nipro.ca | service@nipro.ca

NIPRO CORPORATION

3-26, Senriokashinmachi, Settsu,
Osaka, 566-8510, Japan

© 2023 Nipro Canada. All Rights Reserved.
ParaSil_Rev02
2023-08