




NIPRO | PUSHBAN




Adhesive Pressure Bandage



Before wet swelling

FEATURES



After wet swelling

- The new Nipro PUSHBAN bandage is designed with a special cellulose sponge pad on an adhesive strip.
- The sponge pad is designed to cover the entire wound site, creating a sterile, secured and clean environment.
- This new bandage also expands on contact with blood, which will create additional pressure on the wound site and assist with controlling the bleeding.
- This pressure bandage is ideal for use on hemodialysis access sites, post treatment.
- This pressure bandage is LATEX-FREE

Order Information

Description	Bandage	Size Sponge	Pieces/Box	Box/Case	Re-Order Number
Pushban XL	30x72mm	15x20mm	500	16	PB + XL

Distributed in Canada by:

CardioMed Supplies Inc.

199 Saint David Street
 Lindsay, Ontario K9V 5K7
 Ph: 705-328-2518 Fax: 705-328-9747
 mail@cardiomed.com www.cardiomed.com