

MORE OPTIONS. BETTER OUTCOMES.

The Computed Tomography (CT) compatible polyurethane catheter construction offers radiopacity, durability, and strength.

Simple grip luers are easy on the hands, extremely durable, crack resistant, and won't damage your gloves.

Our super-low profile hub allows for both sutured and sutureless securement.



Easy legible identification on clamps and extensions that indicates power injectability and maximum flow rate.



CardioMed CT Peripherally Inserted Central Catheters (PICC)

Catalogue Number

PICC Description

STANDARD KITS

CMPC-4S60K	4F x 60cm Single Lumen PICC Kit
CMPC-5S60K	5F x 60cm Single Lumen PICC Kit
CMPC-5D60K	5F x 60cm Double Lumen PICC Kit
CMPC-6D60K	6F x 60cm Double Lumen PICC Kit
CMPC-6T60K	6F x 60cm Triple Lumen PICC Kit

MAXIMUM BARRIER TRAYS

CMPC-4S60MBT	4F x 60cm Single Lumen PICC Maximum Barrier Tray
CMPC-5S60MBT	5F x 60cm Single Lumen PICC Maximum Barrier Tray
CMPC-5D60MBT	5F x 60cm Double Lumen PICC Maximum Barrier Tray
CMPC-6D60MBT	6F x 60cm Double Lumen PICC Maximum Barrier Tray
CMPC-6T60MBT	6F x 60cm Triple Lumen PICC Maximum Barrier Tray

All Kits Contain The Following:

CardioMed CT PICC	Safety Scalpel #11
Introducer Needle	Luer Lock Syringe
Catheter Securement Device	Patient Chart Sticker
Patient Information Guide	Tear-Away Introducer
Directions for Use	MicroCLAVE Connectors
70cm Guidewire	24" Measuring Tape

Maximum Barrier Trays Also Include:

Patient Cover	Surgical Gown
Surgical Mask	Bouffant Cap

Non-Latex Surgical Gloves



199 St. David Street
Lindsay, Ontario, Canada
K9V 5K7

1-800-387-9757

Learn more at:

www.CardioMed.com

