



# CardioMed™

Hemoconcentrator Kits

---

## Hemoconcentrator Kits

Nipro Canada's CardioMed™ Hemoconcentrator Kits combine our broad experience in building cardiovascular/perfusion kits with a high-quality Nipro hemoconcentrator membrane.



- ✓ Press-fit Hansen connectors are provided on the kit tubing for making ultrafiltrate and dialysate port connections.
- ✓ Rapid hemoconcentration rates. Virtually eliminates the need for a vacuum source.
- ✓ Biocompatible hemoconcentrator. Housing and fibers are not made with BPA or Bis (2-ethylhexyl) phthalate (DEHP).
- ✓ Hemoconcentrator membrane is made of POLYNEPHRON™, a polyethersulfone (PES) material, made exclusively by Nipro.
- ✓ Kits are available in ultrafiltrate and dialysate-enhanced configurations.

## Ordering Information

Product Code	Description	Membrane Surface Area (m <sup>2</sup> )	Qty per Box
CM9400	Ultrafiltration Tubing Kit	N/A	72
CM-9400-11H	Ultrafiltration Tubing Kit with Hemofilter, 11H	1.1	20
CM9405	Concentrator Tubing Kit	N/A	72
CM-9405-11H	Ultrafilter w/Tubing Kit, 11H	1.1	20
CM-9405-15H	Ultrafilter w/Tubing Kit, 15H	1.5	20
CM-9410	Dialysate Infusion line	N/A	10
CM9504	Hemoconcentrator Tubing Set with 4000mL Collection Bag	N/A	40
CM-9504-09H	Hemoconcentrator Set 0.9m <sup>2</sup> BPA/DEHP Free Elisio Membrane	0.9	20
CM-9504-11H	Hemoconcentrator Set 1.1m <sup>2</sup> BPA/DEHP Free Elisio Membrane	1.1	20
CM-9504-DL	Dialysate Infusion Line / Tubulure Pour L'Infusion De Dialysat	N/A	10



### NIPRO CANADA CORPORATION

199 St. David Street, Lindsay, Ontario, Canada K9V 5K7  
 Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518  
 www.nipro.ca | service@nipro.ca

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE

### NIPRO CORPORATION

3-26, Senriokashinmachi, Settsu,  
 Osaka, 566-8510, Japan

© 2024 Nipro Canada. All Rights Reserved.

HemoconcentratorKits\_Rev03  
 2024-03