

Ultrafil Arterial Filter

Product Data Sheet

| | |
|-------------------------|------------------------|
| Specification | Adult |
| Housing Material | Polycarbonate |
| Filter Fiber Material | Polyester |
| Prime Volume | 179ml |
| Filter Pore Size | 40 µm |
| Blood Ports | 3/8" (10mm) |
| Maximum Blood Flow Rate | 8.0 l/min |
| Filtering Surface Area | Net 500cm ² |
| Max. Operating Time | 6hr |
| Sterilization | EtO Sterilization |

Other Features

- Low Pressure Drop
- Easy to de-bubble
- Very low resistance, no blood trauma produced
- Single Use Device
- Non-Toxic and Pyrogen-free
- Not made with any natural rubber product

Caution

The device must be stored in a dry place with temperature between 0 and 60°C, away from heat sources or bright light.

The device requires proper de-bubbling before use.



CardioMed Supplies Inc.

199 Saint David Street
Lindsay, Ontario, Canada
K9V5K7

www.cardiomed.com

Tel: 1-705-328-2518

Fax: 1-705-328-9747

Toll-free: 1-800-387-9757

mail@cardiomed.com