




***NIPRO*SET™**

B L O O D T U B I N G



**The right line for the right system.
*Every time.***



D
I
A
L
Y
S
I
S
P
R
O
D
U
C
T
S




NIPROSET™


B L O O D T U B I N G



The **NIPRO CONNECT™ PROGRAM** is our commitment to provide every dialysis customer with products of unmatched value through a continuous quality improvement process. In addition, it assures that service, education and the required technical support are provided to each of our business partners.




- **NiproSet™** clear PVC tubing is designed to withstand pressure and manipulation without compromising flow rate accuracy and consistency, regardless of equipment used.
- Laboratory testing demonstrated that when compared with two leading brands of blood tubing, **NiproSet™** delivered an accurate and consistent flow rate over time.



- **NiproSet™** blood lines have demonstrated excellent performance under conditions of higher than normal flow rates and pressures in the dialysis circuit.
- No failures were reported in the blood lines tested in any of the extreme upper limits. Blood and foam never entered the pressure monitor lines and no strike-through occurred.


For the right match, the right flow rate... First time, every time.

All **NiproSet™** tubing is designed to enhance safety and efficiency, with all sets individually pre-labeled based on equipment compatibility. To help ensure the quality you expect, Nipro manufactures a *complete line* of all tubing components, which are tested according to FDA specifications.




MEDICATION PORTS

- Large finger guard helps protect the user from accidental needle sticks; venous line features multiple medication ports for optimum versatility and choice in medication delivery.




COLOR-CODED CLAMPS

- Patient safety is enhanced with clamps that are color-coded for confident identification of the appropriate line.



TRANSDUCER PROTECTORS

- Factory-applied torque helps prevent loosening or disconnecting of transducer; engineered to maintain monitoring while decreasing exposure to the dialysis machine.



DRIP CHAMBERS

- Drip chambers feature a versatile medication port and, for an extra measure of safety, are anti-foaming to minimize splashing.

NIPROSET™

B L O O D T U B I N G

**The right line for the right system.
Every time.**

- All **NiproSet™** tubing features *built-in* safety and efficiency, with all sets pre-labeled based on equipment compatibility and tested according to FDA specifications.
- Laboratory tests demonstrated an accurate and consistent flow rate over time and excellent performance under extreme conditions.
- The **NIPRO CONNECT™ PROGRAM** is our commitment to provide consistently outstanding service, including educational programs that assist in staff training and development.



ORDER SPECIFICATIONS

Reorder Code	Equipment	Pump Segment	Arterial Drip Chamber
A201/V801	Althin 1000/Tina	6.35	Pre-pump
A202/V801	Althin 1000/Tina	6.35	Post-pump
A203/V801	Althin 1000/Tina	8.0	Pre-pump
A204/V801	Althin 1000/Tina	8.0	Post-pump
A219/V803	Fresenius 2008	6.35	Pre-pump
A209/V803	Fresenius 2008	8.0	Pre-pump
A210/V803	Fresenius 2008	8.0	Post-pump
A211/V804	Baxter 550/1550 Meridian	8.0	Post-pump
A212/V804	Baxter Meridian	8.0	Pre-pump
A213/V805	Universal Line	6.35	None
A214/V805	Universal Line	8.0	None
A215/V805	Universal Line	6.35	Pre-pump Pillow
A216/V805	Universal Line	8.0	Pre-pump Pillow
A217/V806	B. Braun	8.0	Post-pump

For additional **NiproSet** configurations, contact us directly at 305-599-7174.

Please consult product labels and inserts for indications, contraindications, hazards, warnings, cautions and directions for use.



NIPRO

NIPRO MEDICAL CORPORATION
3150 N.W. 107 Avenue,
Miami, FL 33172
(305) 599-7174
Fax: (304) 599-8454

For additional information on this product and/or the **NIPRO CONNECT™ PROGRAM**, call us toll-free at 1-888-647-7698/ext. 405.

© 2004, Nipro Medical Corp. Printed in U.S.A.
NiproSet is a trademark of Nipro Medical Corp.
0003 APA 8/04 1M