


Perfusion Accessories



Disposable myocardial needle temperature probe


- ✓ 3.66m (12 ft) long flexible insulated wire with 1/4" right angle plug
- ✓ Compatible with YSI 400 monitors
- ✓ Plugs directly into most temperature monitors
- ✓ Low profile needle housing
- ✓ Fully disposable, eliminating the need for adaptor cables



Non-Adjustable Safety Vent Valve

Based on advanced design, this revolutionary vent valve represents a dramatic improvement in performance, reliability and safety.

- ✓ Clear construction allows complete interior visibility
- ✓ Fits 1/4" tubing



Tubing Organizer

- ✓ Various slots from 7mm to 17mm for different tube sizes.

Perfusion Accessories

Ordering Information

Item	Description	Packaging
2123	Disposable myocardial needle temperature probe (21G x 15 mm)	10 per box
CM-35NAV	One-way vent valve	25 per box
CM-311	Tubing organizer with 10 tube positions	50 per box
CM-311-NS	Tubing organizer with 10 tube positions, non-sterile	100 per box



Quality Management System
Registered to ISO 13485:2003 under CMDCAS

PRODUCT LISTING, DESCRIPTION, AND/OR SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE



CardioMed Supplies Inc.

199 Saint David Street Tel: 705.328.2518
Lindsay, Ontario, Canada Fax: 705.328.9747
K9V5K7 Toll-free: 800.387.9757
www.cardiomed.com mail@cardiomed.com

ACC_PerfusionAccessories_Rev6.0