


Clinical Accessories



MPC200 Tubing Clamps

Available in yellow or blue and made of plastic, CardioMed tubing clamps are appreciated for their durability and cost effectiveness



CM-4516 Cannula Clamps

Fabricated from high grade stainless steel, this "Bulldog" cannula clamp provides consistent clamping pressure, time after time

CM-5015

Recirculating Connector

Made of sturdy plastic, with an attached dust cover, this connector can be used for priming or recirculation



MPC-250 Fistula Clamp

Constructed of durable plastic and metal components, this clamp adjusts to fit the vast majority of patients, generating 1–1.5 lbs of pressure



NIPRO
CANADA

Live Longer. Live Better.

Clinical Accessories

Ordering Information

Code	Description	Packaging
MPC200	Tube occluding forcep (specify colour in blue or yellow)	One hundred (100) per bag
CM-4516	CMS Cannula clamp	Six (6) per bag
CM-5015	CMS Recirculating connector	Fifty (50) per box
MPC-250	Fistula clamp	Ten (10) per bag

Item CM-5015 (recirculating connector) is ETO sterilized

