



DROP-IT™

Quick Release Needle Holder



NIPRO
CANADA

Live Longer. Live Better.

www.nipro.ca

DROP-IT™

Quick Release Needle Holder



Benefits Offered

- Compatible with Vacutainer®¹ or equivalent multi-sample blood collection needles, multi-sample luer adapters, Safety-Lok™¹ wing sets, and blood collection tubes.
- Innovative quick-release mechanism allows for the safe removal of contaminated sharps.
- Can be used up to 100 times.
- Contributes to significant cost savings over single-use disposable blood collection tube holders.
- Transparent holder body facilitates optical veni-puncture control. Holder body made from medical grade material.

Ordering Information

Product Code	Description	Quantity per Box
4550	Drop-It™ Quick Release Needle Holder	50

¹ VACUTAINER® and Safety-Lok™ are trademarks of Becton Dickinson and Company

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE



NIPRO CANADA CORPORATION

199 St. David Street, Lindsay, Ontario, Canada K9V 5K7
Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518
www.nipro.ca | service@nipro.ca

NIPRO CORPORATION

3-26, Senriokashinmachi, Settsu,
Osaka, 566-8510, Japan

© 2024 Nipro Canada. All Rights Reserved.

Drop-It_Rev02
2024-06