

Clinical accessories

Tubing clamps



MPC200

Available in yellow or blue and made of plastic, CardioMed tubing clamps are appreciated for their durability and cost effectiveness

Cannula clamp

CM-4516



Fabricated from high grade stainless steel, this "Bulldog" cannula clamp provides consistent clamping pressure, time after time

Recirculating connector


Made of sturdy plastic, with an attached dust cover, this connector can be used for priming or recirculation



CM-5015

Fistula clamp

MPC-250



Constructed of durable plastic and metal components, this clamp adjusts to fit the vast majority of patients, generating 1–1.5 lbs of pressure

CardioMed
Life-sustaining solutions

Clinical accessories

ORDERING INFORMATION

ITEM	DESCRIPTION	PACKAGING
MPC200	Tube occluding forcep (specify colour in blue or yellow)	One hundred (100) per bag
CM-4516	CMS Cannula clamp	Six (6) per bag
CM-5015	CMS Recirculating connector	One hundred (100) per box
MPC-250	Fistula clamp	Ten (10) per bag
MPC-255	Fistula clamp cover – “Ultra grip”	Twenty-five (25) per bag

Except where otherwise noted, all products are latex-free

Item CM-5015 (recirculating connector) is ETO sterilized

Unless otherwise indicated, products listed are not sterilized



CardioMed Supplies Inc.

199 Saint David Street,
Lindsay, ON K9V 5K7 Canada
Tel: 705.328.2518 Toll-free: 800.387.9757
Fax: 705.328.9747
Web: www.cardiomed.com Email: mail@cardiomed.com