


# Guide wire with heparin coating

## Firm yet flexible

Ideally suited for angioplasty and angiography, CardioMed's guide wires come with a heparin coating that enhances their effectiveness by acting as an anticoagulant.

Flexible enough to be inserted in any blood vessel, CardioMed's guide wires serve to safely direct a catheter during its insertion in a blood vessel.



**CardioMed**  
Life-sustaining solutions

# Guide wire with heparin coating

## ORDERING INFORMATION

ITEM	DESCRIPTION	SPECIFICATIONS
CM-018H12	Guidewire with heparin coating	12.3 cm x .018" with straight end

Except where otherwise noted, all products are latex-free

Unless otherwise indicated, products listed are ETO sterilized



**CardioMed Supplies Inc.**  
199 Saint David Street,  
Lindsay, ON K9V 5K7 Canada  
Tel: 705.328.2518 Toll-free: 800.387.9757  
Fax: 705.328.9747  
Web: [www.cardiomed.com](http://www.cardiomed.com) Email: [mail@cardiomed.com](mailto:mail@cardiomed.com)