



Micro-Solutions for the Neonate

The **NEODRAW** significantly improves arterial line blood sampling from a PAL or UAC catheter to a blood collection device.



Improvements for the Clinician

Ease of Use

- Reduces line setup time
- Split-septum T-connector is the only in-line device needed for arterial line access
- Eliminates need for flushing PAL lines
- Reduces number of line access for blood draw procedure to one
- Standardizes amount of blood waste for every PAL draw versus open drip method

Improves Line Maintenance

- Reduces added bulky components in the bed area used for closed arterial sampling Reduces trouble shooting by elimination of extra in-line components required for closed arterial blood sampling
- No impact on frequency response or wave form from added in-line components

Clinician Safety

- Minimizes blood exposure during blood sampling procedure
- The Hummi Micro-Draw device provides needleless access for blood sampling
- Eliminates open drip method for PAL arterial blood sampling



Ordering Information

Item	Description
CM-ND01	NEODRAW Closed Micro-Draw Blood Transfer Device
Neodraw are packaged 25 per box.	

Improvements for the Neonatal Patient with PAL Line Standardizes the amount of residual blood waste to less than 0.20mL

- Maintains closed system access for drawing blood samples
- Eliminates open drip method of blood for PAL blood sampling
- Eliminates need for flushing PAL line after blood sampling
- PAL blood draw procedure done without positive or negative pressure applied to radial artery
- Fully disposable after single use
- Use only when needed. Eliminates excessive in-line closed blood draw kit costs
- Eliminates set-up time required for use of in-line for closed blood sampling systems

Improvements for the Neonatal Patient with UAC Line

- Maintains closed system access for drawing blood samples
- Reduces number of line access to one per blood draw
- Reduces clearance volume of blood/saline mix to approximately 0.5mL
- Reduces umbilical catheter flush volume after blood draw to approximately 0.5mL Fully disposable after single use
- Use only when needed. Eliminates excessive in-line closed blood draw kit costs
- Eliminates set-up time required for use of in-line for closed blood sampling systems

Technical Specifications

- Accepts luer lock / luer slip devices
- Access with any type of syringe
- Needleless access device



CardioMed Supplies Inc.
199 Saint David Street
Lindsay, Ontario, Canada
K9V 5K7
www. cardiomed.com

Tel: 705.328.2518
Fax: 705.328.9747
Toll-free: 800.387.9757
mail@cardiomed.com