

# CentralCard Vivocath Catheters

## Performance and reliability



CentralCard vivocath radiopaque silicone catheters are engineered for convenience.

Available in single, double and triple lumen configurations in a variety of French sizes.

Designed to support a broad range of treatments cost-effectively.

### Features:

- soft silicone composition that conforms well to vessel anatomy and enhances patient comfort
- compatible with blood products, TPN and a wide range of chemotherapeutic agents in order to perform in a variety of procedures
- integral polyester cuff promotes tissue growth at the exit site, securing it in place and creating a barrier against bacterial infection



Life-sustaining solutions

# CentralCard Vivocath Catheters

## Ordering Information

Item	Description	Lumen	Size
CM-VC64SK	vascular access catheter tray	Single	6.4F x 55cm
CM-VC7SK	vascular access catheter tray	Single	7F x 55cm
CM-VC7DK	vascular access catheter tray	Double	7F x 55cm
CM-VC9SK	vascular access catheter tray	Single	9.5F x 55cm
CM-VC9DK	vascular access catheter tray	Double	9.5F x 55cm
CM-VC12DK	vascular access catheter tray	Double	12.5F x 55cm
CM-VC12TK	vascular access catheter tray	Triple	12.5F x 55cm

CentralCard Vivocath catheters are packaged 5 per box.

CentralCard Vivocath kits include:

- CentralCard Vivocath catheter
- 19G x 2.75" echogenic intro needle
- .035" x 70cm J3 guidewire
- 8" tunneler, soft
- peel-away sheath dilator or peelable introducer
- 1, 2 or 3 needleless swabable adapters

Through continuous innovation, the CardioMed family of CentralCard Vivocath products has grown to include implantable ports, access needles and cap systems.

As with all other CardioMed products, CentralCard Vivocath catheter kits are designed to provide high-quality, cost-effective solutions to real clinical needs and set the standard for patient care.

Non-adhesive catheter repair kits are available for all configurations and make it possible to rapidly restore access for any patient.



**CardioMed Supplies Inc.**  
199 Saint David Street  
Lindsay, Ontario, Canada  
K9V 5K7  
www.cardiomed.com

Tel: 705.328.2518  
Fax: 705.328.9747  
Toll-free: 800.387.9757  
mail@cardiomed.com