

Precise balloon inflation device

Precise for an accurate inflation



Use for the inflation of angioplasty devices requiring dilation.

Features:

- braided high pressure tubing with rotating male luer
- analog luminescent pressure gauge
- comfortable grip handle

Precise balloon inflation device

Ordering Information

Item	Description	Quantity/Box
RS-1011	20ml balloon inflation device	4
RS-1011SK	20ml balloon inflation device with angioplasty accessory components kit (screw-type y connector)	4
RS-1011PK	20ml balloon inflation device with angioplasty accessory components kit (push-type y connector)	4

Angioplasty accessory components kit includes:

- 3-way stopcock
- torque device
- intro needle
- Y connector



Torque device



CardioMed Supplies Inc.
199 Saint David Street
Lindsay, Ontario, Canada
K9V 5K7
www. cardiomed.com

Tel: 705.328.2518
Fax: 705.328.9747
Toll-free: 800.387.9757
mail@cardiomed.com