


Extracorporeal connectors

Top quality connectors for discerning perfusionists

Move your operating room to the best in tubing connectors. CardioMed's exclusive line of extracorporeal connectors are ideal for the discerning perfusionist.









- Enhanced design reduces the incidence of blood hemolysis
- Finger collar ensures a firm grip during placement
- Available with or without a luer lock and offered in a wide range of sizes
- Precision design reduces blood turbulence
- Clear polycarbonate construction allows maximum visibility
- Each connector is supplied with a cap to facilitate sterile technique



CardioMed
Life-sustaining solutions

Extracorporeal connectors

ORDERING INFORMATION

CODE	DESCRIPTION	SIZE	CONFIGURATION
C-656	Y connector without luer lock	1/4" x 1/4" x 1/4"	
C-657	Y connector without luer lock	3/8" x 3/8" x 3/8"	
C-658	Y connector without luer lock	1/2" x 1/2" x 1/2"	
C-659	Y connector without luer lock	3/16" x 3/16" x 3/16"	
C-660	Y connector without luer lock	3/8" x 1/4" x 3/8"	
C-661	Y connector without luer lock	1/2" x 3/8" x 3/8"	
C-662	Y connector without luer lock	3/8" x 1/4" x 1/4"	
C-663	Y connector without luer lock	1/2" x 3/8" x 1/2"	
C-697	Y connector without luer lock	1/2" x 1/4" x 1/2"	
C-664	Y connector with luer lock	1/4" x 1/4" x 1/4"	
C-665	Y connector with luer lock	3/8" x 3/8" x 3/8"	
C-666	Y connector with luer lock	1/2" x 1/2" x 1/2"	
C-667	Y connector with luer lock	3/16" x 3/16" x 3/16"	
C-668	Y connector with luer lock	3/8" x 1/4" x 3/8"	
C-669	Y connector with luer lock	1/2" x 3/8" x 3/8"	
C-670	Y connector with luer lock	3/8" x 1/2" x 1/2"	
C-671	Y connector with luer lock	1/2" x 3/8" x 1/2"	
C-679	Y connector with luer lock	3/8" x 1/4" x 1/4"	
C-672	Straight connector without luer lock	1/4" x 1/4"	
C-673	Straight connector without luer lock	3/8" x 3/8"	
C-674	Straight connector without luer lock	1/2" x 1/2"	
C-675	Straight connector without luer lock	3/16" x 3/16"	
C-684	Straight reducer without luer lock	3/8" x 1/4"	
C-685	Straight reducer without luer lock	1/2" x 3/8"	
C-686	Straight reducer without luer lock	1/4" x 3/16"	
C-680	Straight connector with luer lock	1/4" x 1/4"	
C-681	Straight connector with luer lock	3/8" x 3/8"	
C-682	Straight connector with luer lock	1/2" x 1/2"	
C-683	Straight connector with luer lock	3/16" x 3/16"	
C-676	Straight reducer with luer lock	3/8" x 1/4"	
C-677	Straight reducer with luer lock	1/2" x 3/8"	
C-678	Straight reducer with luer lock	1/4" x 3/16"	
C-689*	1/4" to IV spike connector	1/4"	
C-690	1/4" to male luer lock connector	1/4" with MLL	

* Item C-689 is non-sterile and packaged 25 per box

Unless otherwise indicated, products listed are ETO sterilized and packaged 50 per box

CardioMed Supplies Inc., 199 Saint David Street, Lindsay, ON K9V 5K7 Canada

Tel: 705.328.2518 Toll-free: 800.387.9757 Fax: 705.328.9747 Web: www.cardiomed.com Email: mail@cardiomed.com