


# CentralCard Novacath Venous Catheter Kits



## Safe and reliable

CentralCard NovaCath Venous Catheter is a safe and reliable option for intravenous therapies and monitoring of venous pressure in critical patients.

It is provided as part of a standard kit for catheter insertion via the Seldinger technique.

### Features:

- Catheter body is a biocompatible, radio-opaque, kink resistant PU material.
- Specially designed soft tip prevents endothelial damage during catheter insertion.
- Color-coded extension tubes for quick identification of lumen.
- Incremental markings on the catheter body for quick and easy verification of insertion depth.
- Needleless swabable adapters included with each kit

### Optional

#### Raulerson Bulb Introducer:

- Allows one-step single handed guidewire insertion.
- Reduces the risk of air embolism.
- Prevents blood loss during catheter insertion.
- Minimizes risk of contamination of sterile procedure area.
- Reduces the risk of CRBSI.
- Sold separately.



Life-sustaining solutions

# CentralCard Novacath Venous Catheter Kits

## Ordering Information

Item	Description	Item	Description
CMK-14S15	14G x 15cm single lumen Venous Catheter kit	CMK-4D15	4F x 15cm double lumen Venous Catheter kit
CMK-14S20	14G x 20cm single lumen Venous Catheter kit	CMK-4D20	4F x 20cm double lumen Venous Catheter kit
CMK-14S30	14G x 30cm single lumen Venous Catheter kit	CMK-4D30	4F x 30cm double lumen Venous Catheter kit
CMK-14S45	14G x 45cm single lumen Venous Catheter kit	CMK-5D10	5F x 10cm double lumen Venous Catheter kit
CMK-16S15	16G x 15cm single lumen Venous Catheter kit	CMK-5D13	5F x 13cm double lumen Venous Catheter kit
CMK-16S20	16G x 20cm single lumen Venous Catheter kit	CMK-5D15	5F x 15cm double lumen Venous Catheter kit
CMK-16S30	16G x 30cm single lumen Venous Catheter kit	CMK-5D20	5F x 20cm double lumen Venous Catheter kit
CMK-16S45	16Ga. x 45cm single lumen Venous Catheter kit	CMK-5D30	5F x 30cm double lumen Venous Catheter kit
CMK-18S10	18Ga. x 10cm single lumen Venous Catheter kit	CMK-7D15	7F x 15cm double lumen Venous Catheter kit
CMK-18S15	18Ga. x 15cm single lumen Venous Catheter kit	CMK-7D20	7F x 20cm double lumen Venous Catheter kit
CMK-18S20	18Ga. x 20cm single lumen Venous Catheter kit	CMK-7D30	7F x 30cm double lumen Venous Catheter kit
CMK-18S30	18Ga. x 30cm single lumen Venous Catheter kit	CMK-55T08	5.5F x 8cm triple lumen Venous Catheter kit
CMK-20S10	20Ga. x 10cm single lumen Venous Catheter kit	CMK-55T13	5.5F x 13cm triple lumen Venous Catheter kit
CMK-20S13	20Ga. x 13cm single lumen Venous Catheter kit	CMK-55T15	5.5F x 15cm triple lumen Venous Catheter kit
CMK-20S15	20Ga. x 15cm single lumen Venous Catheter kit	CMK-55T20	5.5F x 20cm triple lumen Venous Catheter kit
CMK-20S20	20Ga. x 20cm single lumen Venous Catheter kit	CMK-55T30	5.5F x 30cm triple lumen Venous Catheter kit
CMK-20S30	20Ga. x 30cm single lumen Venous Catheter kit	CMK-7T15	7F x 15cm triple lumen Venous Catheter kit
CMK-24S8	24Ga. x 8cm single lumen Venous Catheter kit	CMK-7T20	7F x 20cm triple lumen Venous Catheter kit
CMK-24S13	24Ga. x 13cm single lumen Venous Catheter kit	CMK-7T30	7F x 30cm triple lumen Venous Catheter kit
CMK-24S15	24Ga. x 15cm single lumen Venous Catheter kit	CMK-85Q15	8.5F x 15cm 8.5F quad lumen Venous Catheter kit
CMK-24S20	24Ga. x 20cm single lumen Venous Catheter kit	CMK-85Q20	8.5F x 20cm 8.5F quad lumen Venous Catheter kit
CMK-4D05	4F x 5cm double lumen Venous Catheter kit	CMK-85Q30	8.5F x 30cm 8.5F quad lumen Venous Catheter kit
CMK-4D08	4F x 8cm double lumen Venous Catheter kit	17-6001	18Ga. x 7cm Raulerson bulb introducer set
CMK-4D13	4F x 13cm double lumen Venous Catheter kit	CentralCard Novacath Venous Catheter kits are packaged 10 per box.	

CentralCard Novacath Venous Catheter kits include:

- CentralCard Novacath venous catheter
- Vessel dilator
- Guidewire
- Introducer needle
- Movable fixation wing
- #11 long handle scalpel
- Needleless swabable adapters
- 5cc syringe with MLL
- 22G needle
- 25G needle
- 5cc Raulerson syringe\*\*

\*\* (not included in 4F double lumen or 20G/24G kits)



**CardioMed Supplies Inc.**  
 199 Saint David Street  
 Lindsay, Ontario, Canada  
 K9V 5K7  
[www.cardiomed.com](http://www.cardiomed.com)

Tel: 705.328.2518  
 Fax: 705.328.9747  
 Toll-free: 800.387.9757  
[mail@cardiomed.com](mailto:mail@cardiomed.com)