

INTRODUCING

# EXOFLEX™



## The manifold that morphs with your clinical demands

ExoFlex™ with Q2® Technology is a new manifold designed to span departments and accommodate the need of being **rigid or flexible - in seconds.**

- Adaptable between departments and clinicians
- Reduce risk of HA-BSI with the Q2 swabable, split septum valve<sup>1</sup>
- Low residual volume/dead space
- Non-hemolytic Q2 swabable valve for repeated access<sup>2</sup>
- Integrated check valves ensure one-way infusion to the patient; preventing retrograde flow
- High flow rates for bolus needs
- Compatible with ISO 594 / 80369-7 disinfectant caps containing 70% IPA<sup>2</sup>
- ExoFlex decreases line manipulation reducing associated errors and infection risks



**ExoFlex. Spanning departments without breaking the line.**

Contact us for a sample or for more details.

Exclusively Distributed by:



199 St. David Street | Lindsay, Ontario, Canada | K9V-5K7  
705-328-2518 | 800-387-9757  
www.cardiomed.com



 **QUEST** Medical, Inc.

custserv@questmedical.com | +1.800.627.0226  
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com

1 CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services; 2 Data on file  
© 2019 Quest Medical, Inc. 2020-28

INTRODUCING

# EXOFLEX™



## The manifold that morphs with your clinician's demands

ExoFlex™ with Q2® Technology is a new manifold designed to span departments and accommodate the need of being **rigid or flexible - in seconds.**

- Assists in waste reduction initiatives
- Reduce risk of HA-BSI with the Q2 swabable, split septum valve<sup>1</sup>
- Reduce soft costs as the product is assembled, sterilized and ready for use
- Reduce SKUs
- ExoFlex decreases line manipulation reducing associated errors and infection risks



**ExoFlex. Spanning departments without breaking the line.**

Contact us for a sample or for more details.

Exclusively Distributed by:



199 St. David Street | Lindsay, Ontario, Canada | K9V-5K7  
705-328-2518 | 800-387-9757  
www.cardiomed.com

 **QUEST** Medical, Inc.

custserv@questmedical.com | +1.800.627.0226  
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com



<sup>1</sup> CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services; © 2019 Quest Medical, Inc. 2020-28