



Optimum Flow Rates

Maximize Delivery

High Flow Stopcock

9 French Internal Diameter
Sterile STV01 | Non-Sterile STV01BNS

Priming Volume: 0.5 mL
Flow Rate: 1,360 mL/min

Capless High Flow Stopcock prevents bottle-necking of fluid and drugs being delivered, allowing for maximized delivery of infusates.



High flow rate¹



Minimize infection risk²



Alcohol and lipid resistant¹



High Flow Stopcock

Maximize Delivery of Drugs and Fluids

Quest's Capless High Flow Stopcock

- 9 French inner diameter prevents fluids from bottle-necking and allows high flow rates for maximum delivery
- Whole solution delivery ensured via Q2[®] Valves with no internal componentry
- Minimize infection risk with closed system
- Alcohol-resistant housing for peace of mind use



Quest Split Septum
Needleless Connectors

