








SRV™ DATA SHEET

-  **Single Use**
-  **Connects using standard Luer Connectors**
-  **Makes return to Treatment fast, simple, & clean**
-  **Proprietary Flow Path Seals when Valve is**
-  **Activated under predetermined tension**



SAFETY RELEASE VALVE™ - Performance Specifications

Luer Connectivity	Yes
Single Use Only	Yes
Materials	Silicone, Polycarbonate
Priming Volume	0.182 mL
Residual Volume following Activation	0.058 mL
Non-Metallic	Yes
Flow Path	In-line proprietary design
Latex/DEHP Free	Yes
Disconnect Tension	Max 3.25 lbs

SAFETY RELEASE VALVE™ - Order Information

Reference Ordering Number	1812-B300 Rev D
Packaging	Individual pkgd unit
Sold As Case	100 units