



# rescuCLAMP<sup>®</sup>

Catheter Repair Clamp

---



[www.nipro.ca](http://www.nipro.ca)

# rescuCLAMP®

## Catheter Repair Clamp

Nipro's Catheter Repair Clamp allows the clinician to quickly replace a malfunctioning or broken pinch clamp by inserting the tubing through the specially designed side of the repair clamp.



Unique grooves allow the repair clamp to be mounted easily and securely

Guide pin keeps the clamp centered ensuring that the extension tube is pinched evenly



## Ordering information

Product code	Description	Size	Packaging
CM-604	rescuCLAMP®	Fits .118" - .177" OD Tubing	25 per Box
Non Sterile		No Latex	

# rescuCLAMP®

Catheter Repair Clamp



YOUR LINK BETWEEN PATIENT AND CARE



**NIPRO CANADA CORPORATION**

199 Saint David Street, Lindsay, Ontario, Canada K9V-5K7

Toll-Free: 800-387-9757 | Local: 705-328-2518

[www.nipro.ca](http://www.nipro.ca) | [service@nipro.ca](mailto:service@nipro.ca)

**NIPRO CORPORATION**

3-9-3 Honjo-Nishi

Kita-Ku, 27, JP, 531-8510

© 2022 Nipro Canada. All Rights Reserved.  
rescuCLAMP\_Rev03