

CARDIOVASCULAR

INFLATION DEVICES


dolphin



dolphin

The inflation device matching the majority of balloon catheter

- 30 atm / 3040 kPa pressure gauge
- 30 cm³ syringe
- Ergonomic handle



FEATURES & BENEFITS

PRESSURE GAUGE

- 30 atm / 3040 kPa
- suitable for most of the procedures

BODY AND BUTTON

- 30 cm³ syringe
- small to big balloons (PTA & PTCA)
- Square shaped
- stable
- Automatic cam locking system
- easy handling
- Clear material
- easy air bubbles detection and removal

HANDLE

- T-shape handle
- high torque
- easy deflation

REFERENCES

References	Description	Accessories	Packaging
0185ND	30 atm / 3040 kPa inflation device		6/box
0185NA	30 atm / 3040 kPa inflation device	Trefoil 3-way stopcock	
0185PD	30 atm / 3040 kPa inflation device	Y connector push-pull Metallic guidewire introducer Alligator torque device Trefoil 3-way stopcock	
0185QL	30 atm / 3040 kPa inflation device + PTCA set 3 items	Trefoil 3-way stopcock myshell lite Y connector 9F / 3mm, 20cm PVC connection line - 3-way stopcock Metallic guidewire introducer Alligator torque device	

Please contact VYGON for non sterile and OEM products.

Please read carefully the Instructions For Use. The CE conformity assessment has been conducted by GMED (0459).



05/2023 - DB BU CS 18 439 E

FOR FURTHER INFORMATION, PLEASE CONTACT: questions@vygon.com

The specification shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

VYGON – 5, RUE ADELINÉ - 95440 ECOUEN - FRANCE

RECEPTION: +33 (0)1 39 92 63 63 – SERVICE CLIENTS FRANCE: +33 (0)1 39 92 63 81

EXPORT CUSTOMER SERVICE: +33 (0)1 39 92 64 15

www.vygon.com

