



SmartTap[®]

Paracentesis Drainage Trays

Unparalleled safety, exceptional performance

Paracentesis Drainage Trays



Nipro Canada’s sterilized SmartTap® paracentesis trays can be used for diagnostic, therapeutic, or palliative fluid withdrawal. Additional sterilized drainage bags can be purchased standalone.

Benefits

- Dedicated trays provide clinicians with all the components required for safe and effective drainage procedures.
- Closed-system design reduces treatment time and the risk of cross-contamination.
- Both tray kits come with safety syringes with retractable needle.
- 14G catheter with retractable safety needle available (CM-002P-S only).

Paracentesis trays include:

- 4L drainage bag
- 14G catheter over needle with safety (CM-002P-S only)
- 16G needle with 14G over needle catheter with side holes (CM-0020T and CM-8016 only)



- 10cc luer lock syringe
- 5cc safety syringe with 22G x 1 1/2” retractable needle
- 3cc safety syringe with 25G x 1” retractable needle
- Antiseptic Swabstick with 2% w/v Chlorhexidine Gluconate and 70% v/v Isopropyl Alcohol
- 4-way stopcock
- Gauze pads
- Fenestrated drape and underpad

Ordering Information

Product Code	Description	Qty per Box
CM-002P-S	SmartTap® Paracentesis Tray w/3.3cm Catheter, Safety Needle	8
CM-0020T	SmartTap® Paracentesis Drainage Tray w/9cm Catheter	8
CM-0020	SmartTap® Paracentesis Drainage Bag, 4L	25
CM-0022	SmartTap® Paracentesis Drainage Bag, 2L	25
CM-8016	14G x 9cm Teflon Catheter Over Needle w/Side Holes	10

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE

