



VanishPoint® II

Same Retraction and Safety!

SYRINGE IMPROVED FOR YOU!

- Enhanced aspiration and bubble removal
- Improved plunger handle control
- Increased user comfort during administration
- Reduced pressure required to depress plunger

New VanishPoint II plunger handle design provides **increased surface contact for your thumb, improving comfort and control** during aspiration and injection!



- GOOD NEWS!** The new VanishPoint II syringe is being brought to you as the result of customer clinical input. While we are recognized as the world leader in retractable needle technology, we continue to work diligently to improve our innovative products even further.
- We listen carefully to product suggestions from our customers, and the new VanishPoint II syringe is an important result. This new syringe incorporates more efficient aspiration and bubble removal. In addition, we have modified the syringe's plunger handle to a concave design that makes it more comfortable to depress. The advantages of these improvements became evident during

the COVID-19 pandemic. RTI supplied hundreds of millions of VanishPoint syringes that were used to vaccinate American residents against the COVID corona virus. During the pandemic, healthcare workers, pharmacists, and others were performing unusually high numbers of injections per day. As a result of their input, we modified the plunger handle ergonomically so that it would require less thumb pressure to activate the retraction mechanism.

- We proudly introduce to you the VanishPoint II syringe, a product that is improved. Our customers asked for it, we listened, and we responded. **THANKS FOR YOUR INPUT!**



www.retractable.com

rtisales@retractable.com

511 Lobo Lane • Little Elm • Texas • 75068-5295 • USA
Tel: (972) 294-1770 • Fax: (972) 294-4400
Toll Free: 1-888-703-1010



RETRACTABLE
TECHNOLOGIES, INC.