



CARDIOPLEGIA ADAPTERS

CARDIOPLEGIA ADAPTERS

Facilitate optimal heart protection during cardiac surgery with versatile cardioplegia adapters.


These adapters offer unmatched customization for your cardioplegia circuit, enabling you to:

- tailor delivery methods: Switch between antegrade and retrograde cardioplegia with ease.
- accommodate various techniques: Multiple configurations ensure a perfect fit for your surgical approach.
- enhance efficiency: Simplify circuit setup and optimize blood flow management.

Available in a range of configurations, including:

- Vent adapters for air removal
- Multi-perfusion adapters for reaching multiple targets

All Cardioplegia Adapters include standard luer lock connectors.



Ordering Information

Item	Description	Units per Box
SPBC-R	Cardioplegia Y-Adapter w/One Venting Arm	20
SPBC-BR	Cardioplegia Y-Adapter w/One Extended Arm	20
SPBC-41	Cardioplegia Perfusion Adapter Triplicate	20
SPBC-SH	Cardioplegia Perfusion Adapter - Quad	20
SPBC-41BR	Cardioplegia Perfusion Adapter Quad w/Vent Line	20



NIPRO CANADA CORPORATION
 199 St. David Street, Lindsay, Ontario, Canada K9V 5K7
 Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518
 www.nipro.ca | service@nipro.ca

NIPRO CORPORATION
 3-26, Senriokashinmachi, Settsu,
 Osaka, 566-8510, Japan

© 2024 Nipro Canada. All Rights Reserved.
 Cardioplegia Adapters Rev_01 2024-05