



Providing a New Alternative to Hemodialysis Patient Care

Uniquely designed for use in hemodialysis, the **Lyka™ by Quest Medical** brings engineering excellence to patient care, safety, and comfort.

The needleless connector utilizes the Quest split septum with low dead space and is fully flushable to reduce infection risk. With an open fluid path, it provides Ultra High flow rates ideal for your patient's needs. The closed nature of the valve offers protection to the patient against catheter-related blood stream infections (CRBSIs) and helps reduce clinical costs related to CRBSIs.

Lyka Key Features:

- High flow rate of up to 600 mL/min²
- Straight fluid path, low dead space, and fully flushable
- Needleless connector creates a closed-type system to help protect the patient against CRBSIs
- Translucent housing for visual check
- Provides an effective microbiological barrier³
- Approved for 7 day use

Provide a safe and effective solution for patient care during hemodialysis treatments with the new **Lyka by Quest Medical**.¹ To learn more about the newest option in dialysis care, email us at lyka@questmedical.com or call us at **1.800.627.0226**.



Straight fluid path for high flow rates



Allen, Texas | +1.972.390.9800 | +1.800.627.0226
lyka@questmedical.com | www.questmedical.com

Rx ONLY



For over 40 years, Quest Medical, Inc. has been dedicated to providing the highest quality products for patient care. Our products are designed, engineered, and manufactured in the U.S.A.

¹ Not available for sale in the U.S.A.
² Data on file
³ Data on file