

**For Fluid Control and Solution Administration**




17-469N

# STOPCOCKS

# STOPCOCKS

Nipro Canada's Stopcock product family offers the flexibility and reliability expected from health care providers, coming in various pressure, flow, directional, and material configurations.



17-130 shown above

- ✓ The 4-way 17-469 stopcock comes with a polysulfone body capable of being used for up to 96 hours with most common drug solutions such as intralipids.
- ✓ Available in standard, medium, and high-pressure rated configurations.
- ✓ Large bore models are available for greater flow rate.

## Ordering Information

Product Code	Description	Quantity per Box
17-300	One Way Stopcock w/Rotating Male L.L.	50
17-369C	3-Way Stopcock with Caps	50
17-369H	3 Way High Pressure Stopcock - 1200 psi	50
17-469	Lipid Compatible Large-Bore Stopcock	50
17-469N	4-Way, Large-Bore Stopcock w/Swabable Port	100
17-305	1-Way High Pressure Stopcock, White Handle w/Swivel M.L.L. - 1200 PSI	50
17-308	3-Way Medium Pressure Stopcock, White Handle w/Swivel M.L.L. - 500 PSI	50
17-309	1-Way High Pressure Stopcock, White Handle w/Swivel High Pressure Adapter M.L.L. - 1200 PSI	50
17-130	Anesthesia Manifold	25
17-131	Perfusion Manifold	25



### NIPRO CANADA CORPORATION

199 St. David Street, Lindsay, Ontario, Canada K9V 5K7  
Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518  
www.nipro.ca | service@nipro.ca

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE

### NIPRO CORPORATION

3-26, Senriokashinmachi, Settsu,  
Osaka, 566-8510, Japan

© 2025 Nipro Canada. All Rights Reserved.  
Stopcocks\_Rev00  
2025-09