




HEMOCONCENTRATION KITS

HEMOCONCENTRATION KITS

Nipro Canada's Hemoconcentration Kits combine our broad experience in building cardiovascular/perfusion kits with a high-quality hemoconcentrator membrane.



CM-9405-11H shown

- ✓ Rapid hemoconcentration rates. Virtually eliminates the need for a vacuum source.
- ✓ Biocompatible Nipro Elisio™ hemoconcentrator. Housing and fibers are not made with BPA or Bis (2-ethylhexyl) phthalate (DEHP).
- ✓ Press-fit Hansen connectors are utilized for making ultrafiltrate and dialysate port connections.
- ✓ Dialysate-enhanced infusion lines are available as an option.

Ordering Information

Product Code	Description	Quantity per Box
CM-9400	16" Blood Lines - Tubing Kit Only	5
CM-9400-09H	0.9m ² Hemoconcentrator with CM-9400 Tubing Kit	20
CM-9400-11H	1.1m ² Hemoconcentrator with CM-9400 Tubing Kit	20
CM-9400-15H	1.5m ² Hemoconcentrator with CM-9400 Tubing Kit	20
CM-9405	36" Blood Lines - Tubing Kit Only	5
CM-9405-09H	0.9m ² Hemoconcentrator with CM-9405 Tubing Kit	20
CM-9405-11H	1.1m ² Hemoconcentrator with CM-9405 Tubing Kit	20
CM-9405-15H	1.5m ² Hemoconcentrator with CM-9405 Tubing Kit	20
CM-9410	Dialysate Infusion Line for CM-9400 and CM-9405	10
CM-9504	36" Blood Lines - Tubing Kit w/4L Collection Bag Only	5
CM-9504-09H	0.9m ² Hemoconcentrator with CM-9504 Tubing Kit	20
CM-9504-11H	1.1m ² Hemoconcentrator with CM-9504 Tubing Kit	20
CM-9504-15H	1.5m ² Hemoconcentrator with CM-9504 Tubing Kit	20
CM-9504-DL	Dialysate Infusion Line for CM-9504	10
CM-9515	15" Blood Lines - Tubing Kit w/4L Collection Bag Only	5
CM-9515-09H	0.9m ² Hemoconcentrator with CM-9515 Tubing Kit	20
CM-9515-11H	1.1m ² Hemoconcentrator with CM-9515 Tubing Kit	20
CM-9515-15H	1.5m ² Hemoconcentrator with CM-9515 Tubing Kit	20



NIPRO CANADA CORPORATION

199 St. David Street, Lindsay, Ontario, Canada K9V 5K7
 Toll-Free: 1-800-387-9757 | Local: 1-705-328-2518
www.nipro.ca | service@nipro.ca

PRODUCT LISTING, DESCRIPTION AND SPECIFICATION ARE SUBJECT TO CHANGE WITHOUT NOTICE

NIPRO CORPORATION

3-26, Senriokashinmachi, Settsu,
 Osaka, 566-8510, Japan

© 2026 Nipro Canada. All Rights Reserved.
 HemoconcentrationKits_Rev05
 2026-02