



rescuCLAMP®

Catheter Repair Clamp



www.nipro.ca

rescuCLAMP®

Catheter Repair Clamp

Nipro's Catheter Repair Clamp allows the clinician to quickly replace a malfunctioning or broken pinch clamp by inserting the tubing through the specially designed side of the repair clamp.



Unique grooves allow the repair clamp to be mounted easily and securely

Guide pin keeps the clamp centered ensuring that the extension tube is pinched evenly



Ordering information

Product code	Description	Size	Packaging
CM-604	rescuCLAMP®	Fits .118" - .177" OD Tubing	25 per Box
Non Sterile		No Latex	

rescuCLAMP®

Catheter Repair Clamp



NIPRO
CANADA

Live Longer. Live Better.



NIPRO
CANADA

NIPRO CANADA CORPORATION

199 Saint David Street, Lindsay, Ontario, Canada K9V-5K7

Toll-Free: 800-387-9757 | Local: 705-328-2518

www.nipro.ca | service@nipro.ca

NIPRO CORPORATION

3-9-3 Honjo-Nishi

Kita-Ku, 27, JP, 531-8510

© 2022 Nipro Canada. All Rights Reserved.
rescuCLAMP_Rev03