trio-ct™

TRIPLE LUMEN CATHETER

ACUTE DIALYSIS CARE

Ordering Information

TRIO-CT™ BASIC SET (CANADA)



CATHETER SIZE	INSERTION LENGTH	CATHETER TYPE	BASIC SET (5/BOX)
13.5F	12cm	Straight	AS01135121CA
13.5F	15cm	Straight	AS01135151CA
13.5F	20cm	Straight	AS01135201CA
13.5F	24cm	Straight	AS01135241CA
13.5F	30cm	Straight	AS01135301CA

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (2) Dilators (1) Scalpel (1) Adhesive Wound Dressing (1) Attachable Suture Wing (3) Injection Caps (1) Patient ID Card (1) Chart Sticker

Medical Components, Inc.
1499 Delp Dr., Harleysville, PA 19438, U.S.A.
P: 215.256.4201 **F**: 215.256.1787
www.**medcompnet**.com

PN2643 Rev. A 10/18



***Three** lumens...*
***One** catheter...*
Perfect Harmony


medCOMP®

Indicated for hemodialysis and apheresis.
Third lumen provides additional access for multiple therapies.

trio-ct™

TRIPLE LUMEN CATHETER

INTERNAL LUMEN



17 ga lumen is independent from dialysis lumens and exits distal to the dialysis lumens.

CATHETER TIP



Tapered catheter tip design.

KINK RESISTANT



13.5F oval catheter designed for kink resistance.



Power Injectable Lumen

Extension is labeled "Power Injectable."

Third lumen decreases the need for additional needle sticks.


Provides access for infusion of medication, blood draws, central venous pressure monitoring, and power injection of contrast media.

Allows injection of contrast media at a maximum of 5 mL/sec.

Average recirculation rate of 2% for normal flow, and when arterial and venous lumens are reversed

Flow rates of 400 ml/min and a pressure of 80 mmHg for straight 24 cm catheter length

MRI safe



ID tags indicate maximum injection rates

Priming volume shown on catheter clamps

Catheter lengths up to 30 cm

Depth markings aid with catheter insertion to desired depth

medCOMP®

www.medcompnet.com